Platelet scale set

REF 1479

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 0SD, United Kingdom

Tel: +44 (0)191 482 8440

Fax: +44 (0)191 482 8442

Email: info@helena-biosciences.com

Web: www.helena-biosciences.com

HL-2-1568P 2015/10 (5)

Platelet scale set Instructions for use

en

INTENDED PURPOSE

The Platelet scale set kit is intended for use as a calibration material.

Helena Biosciences Europe Platelet scale set is intended to be used to optically calibrate the AggRAM at 650 nm.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *In Vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

COMPOSITION

Component	Content	Description	Preparation
Scale set 1	1 x 7 mL	Scale set 1 is an aqueous solution containing buffers and stabilizers.	Allow solutions to warm to room temperature.
Scale set 2	1 x 3 mL	Scale set 2 additionally contains cross-linked styrene divinylbenzene copolymer microbeads.	Allow solutions to warm to room temperature.

Each kit contains instructions for use.

Each kit contains a lot specific reference values insert.

ITEMS REQUIRED BUT NOT PROVIDED

REF 1487 / 1487PC AggRAM

STORAGE, SHELF-LIFE AND STABILITY

The solutions should be stored at $2^{\circ}\text{--}8^{\circ}\text{C}$, and are stable until the expiration date indicated on the package. DO NOT FREEZE!

Component Signs of Deterioration

Scale set 1 cloudiness may indicate contamination.

Scale set 2 cloudiness may indicate contamination. The Scale set 2 solution should be a whitish, turbid suspension. Any non-white colour may indicate contamination. Avoid freezing, which can cause the microbeads to disintegrate.

PROCEDURE

Method AggRAM (REF 1487 / 1487PC)

Perform Scale set monthly.

1. Remove Scale set 1 and Scale set 2 from refrigerator and allow to come to room temperature ($15^{\circ}\text{--}30^{\circ}\text{C}$).
2. Allow the instrument to warm up for 20 minutes or until the optics reach 37°C .
3. Perform the daily Optical Calibration Check, see section 10.1.1 of the Instrument Operators' Manual.
4. Prepare the instrument to perform the Scale set, see section 10.1.5 of the Operators' Manual and enter the Scale set value provided at the end of this procedure.
5. Place two covers in the cuvette rack.
6. Remove the scale of the filled cuvette. Return the cuvette to the rack.
7. Using the same method, place at least 400 μL of Scale set 2 into the second cuvette.
8. Scale set 1 is checked first, see section 10.1.5 of the Operators' Manual for instructions.
9. Scale set 2 is checked second, see section 10.1.5 of the Operators' Manual for instructions.
10. Values for Scale set 1, and Result (V/O.D.) display and are colour coded to indicate acceptable/unacceptable values. See section 10.1.5 of the Operators' Manual for information on the colour meanings.
11. If the values are unacceptable, repeat the Scale set procedure and if any values remain unacceptable see section 10.1.5 of the Operators' Manual.
12. Print a hard copy of the Scale set values by selecting Print.
13. The Scale set values are automatically saved when the window is exited.

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