REF 5562
REF 5558
REF 5559

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom
Tel: +44 (0)191 482 8440
Fax: +44 (0)191 482 8442
Email: info@helena-biosciences.com
Web: www.helena-biosciences.com

HL-2-1787P 2016/01 (6)

APTT Si L Minus Instructions for use

INTENDED PURPOSE

The APTT Si L Minus kit is intended for carrying out clot based haemostasis assays.

For use in the determination of activated partial thromboplastin times (aPTT), and related coagulation procedures using phospholipid extract and a near-colloidal particle activator. The test system can be used on manual, semi-automated and automated methods. From its origins through the work of Langdell and coworkers¹, and later modified by Proctor and Rapaport², the aPTT is used to detect disorders in the intrinsic coagulation system, which involves coagulation factors VIII, IX, XI, XII, prekallikrein, and high molecular weight kininogen. The aPTT is also used in assays which quantitate these factors and is routinely used for presurgical screening and monitoring of heparin therapy³. Commercially available reagents typically use one of three activators: kaolin, silica, or ellagic acid. In the basic screening test, the aPTT indirectly measures the formation of thrombin by its action on fibrinogen forming the fibrin clot. In the test, citrated test plasma is mixed with aPTT reagent for a specified period of time (typically 5 minutes) at 37°C followed by the addition of pre-warmed (37°C) calcium chloride (0.025 M). Timing is begun from the time of addition of calcium chloride. The time required for clot formation is the aPTT. Clot detection can be by mechanical, manual (tilt tube), or photo-optical measurement.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

COMPOSITION

Component	Content	Description	Preparation
APTT Si L Minus	5 x 5 mL (REF 5562 / REF 5562SLQ) 5 x 10 mL (REF 5558 / REF 5558SLQ) 10 x 10 mL (REF 5559 / REF 5559SLQ)	Reagent contains a near colloidal particle activator (magnesium-aluminium-silicate) for optimum sensitivity to factor deficiencies and to heparin. The reagent also contains phospholipids with buffer and stabilisers.	Bring to room temperature prior to use. Mix well by swirling or inversion prior to use.
Calcium Chloride: 0.025M	5 x 5 mL (REF 5562 / REF 5562SLQ) 5 x 10 mL (REF 5558 / REF 5558SLQ) 10 x 10 mL (REF 5559 / REF 5559SLQ)	The reagent is a 0.025 M solution of calcium chloride.	The reagent is ready for use as packaged.

ITEMS REQUIRED BUT NOT PROVIDED

Any high quality electro-mechanical or photo-optical coagulation instrument designed for performing activated partial thromboplastin times may be used.

STORAGE, SHELF-LIFE AND STABILITY

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label. Store at 2° – 8°C. DO NOT FREEZE. Stable for 30 days after opening. Avoid prolonged heating.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at 2° – 8°C or 18 – 24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at 37°C prior to testing. Do not keep at 37°C for more than 5 minutes. This will minimise the neutralisation of the lupus inhibitor⁴. Erroneous results may be caused by contamination with tissue fluids or stasis. Avoid agitation, air bubbles or foaming. For the effects of commonly administered drugs, refer to Young, *et al.*⁵.

PROCEDURE

Manual Method

- Prewarm well mixed APTT Si L Minus and 0.025 M Calcium Chloride to 37°C.
- Prewarm 0.1 mL test plasma in duplicate to 37°C for 2 minutes.
- Forcibly add 0.1 mL prewarmed APTT Si L Minus to plasma and start timer. Incubate for exactly 5 minutes at 37°C.
- Add 0.1 mL prewarmed 0.025 M Calcium Chloride.
- Note time for clot formation. Report result as aPTT time (seconds).

Automated Method

Refer to the appropriate instrument operator manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application notes.

INTERPRETATION OF RESULTS

The results of the aPTT test should be reported to the nearest 1/10 of a second. The normal range (usually X ± 2 standard deviations) for each individual laboratory should be established. Results greater than the upper limits of the normal range should be considered abnormal and follow-up testing should be performed. Any aPTT values less than the lower limits of the normal range should be repeated on a new blood sample. Short aPTT values may be seen in association with *in vivo* thrombosis (e.g. deep vein thrombosis and disseminated intravascular coagulation).

Heparin monitoring

When monitoring heparin therapy, it is important to construct an *in vitro* reference curve which reflects the average heparin response, since individual patients respond differently to heparin. In general, one can consider the therapeutic range for heparin to be 0.2 to 0.5 units/mL^{6,8}.

The following precautions should be considered when monitoring heparin therapy:

- Time of collection is important, since heparin has an *in vivo* half-life of only 1.5 hours⁷.
- Release of platelet factor 4 (heparin neutralising factor) caused by platelet aggregation or damage during collection, should be avoided. Careful blood collection, proper centrifugation and prompt removal of the platelet poor plasma from the cells will help minimise the release of platelet factor 4.
- Baseline data on each patient's aPTT should be established before therapy, to determine the respective patient aPTT as it relates to the normal range established by the laboratory.
- Different clot detection systems (mechanical, photo-optical, etc.) show variable sensitivities to heparin. The same test system should be used when monitoring heparinised patients.
- Heparin response curves should be reestablished when lot numbers of reagent change and at periodic intervals with the same lot number.
- The curve should also be constructed using the same heparin employed in therapy, to eliminate variables connected with heparins from different sources (e.g. porcine mucosa or bovine lung).

LIMITATIONS

Expected values for the aPTT test will vary from one laboratory to another, depending on the technique used. The method of clot detection, temperature, pH, collection technique, type of anticoagulant and time and method of specimen storage are all very important. Plasma sample collection and storage conditions should be standardised and carefully controlled. Unexpected results should be confirmed by additional tests. Platelet fragments present in a specimen may cause the release of phospholipids, and thus the neutralisation of any lupus inhibitor present in the specimen. The use of specimens with small plasma volumes should be avoided due to possible physiological pH changes. Testing could be affected by several drugs⁵. An increase in the aPTT results may be caused by the administration of diphenhydantoin, heparin, warfarin and radiographic agents⁸. Decreased aPTT values may be seen during the use of oral contraceptives, or male estrogen therapy^{9,10}. Thus, laboratories should establish their own expected values for patients and well defined performance standards for the control.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. Helena Biosciences Europe supplies the following controls available for use with this product:

REF 5188	Routine Control N
REF 5187	Routine Control A
REF 5183	Routine Control SA
REF 5189	Heparin Control
REF 5301	Speciality Assayed Control N
REF 5302	Speciality Assayed Control A

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives using a photo-optical coagulation instrument. Each laboratory should establish its own performance data.

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BIBLIOGRAPHY

- Langdell R, Wagner R, Brinkhous K (1953) Effect of Antihemophilic Factor on One Stage Clotting Tests, *J. Lab. Clin. Med.* **41**: 637
- Proctor R, Rapaport S (1961) The Partial Thromboplastin Time with Kaolin, *Am. J. Clin. Path.* **36**: 212
- Brandt JT and Triplett DA (1981) Laboratory Monitoring of Heparin. Effect of Reagents and Instruments on the Activated Partial Thromboplastin Time, *AJCP*, **76**: 530-537
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Young DS *et al.* Effects of Drugs on Clinical Laboratory Tests, 3rd ed., AAC Press, Washington, D.C., 1990
- Kelton JG and Hirsh J (1980) Bleeding Associated with Antithrombotic Therapy, *Semin Hematol*, **17**: 259-91
- Estes J and Poulin P (1974) Pharmacokinetics of Heparin, *Thromb Diath Haemorr*, **33**: 26-35
- Solomon G, Hilgartner M, Kutt H (1972) Coagulation defects caused by diphenhydantoin, *Neurology*, **22**: 1165-71
- Ambrus JL, Schimert G, Lajos TZ, Ambrus CM, Mink IB, Lassarman HB, Moore RH, Meizer J (1971) Effect of antifibrinolytic agents and estrogens on blood loss and coagulation factors during open heart surgery, *J Med*, **2**: 65-81
- Crowell EB Jr, Clatanoff DV, Kiekhof W (1971) The effect of oral contraceptives on factor VIII levels, *Lab. Clin. Med.* **77**: 551-7

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BIBLIOGRAFIA
<ol style="list-style-type: none">Langdell R, Wagner R, Brinkhous K (1953) Effect of Antihemophilic Factor on One Stage Clotting Tests, <i>J. Lab. Clin. Med.</i> 41: 637 Proctor R, Rapaport S (1961) The Partial Thromboplastin Time with Kaolin, <i>Am. J. Clin. Path.</i> 36: 212 Brandt JT and Triplett DA (1981) Laboratory Monitoring of Heparin. Effect of Reagents and Instruments on the Activated Partial Thromboplastin Time, <i>AJCP</i>. 76: 530-537 Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5 Young DS <i>et al.</i> Effects of Drugs on Clinical Laboratory Tests, 3rd ed., AAC Press, Washington, D.C., 1990 Kelton JG and Hirsh J (1980) Bleeding Associated with Antithrombotic Therapy, <i>Semin Hematol</i>, 17: 259-91 Estes J and Poulin P (1974) Pharmacokinetics of Heparin, <i>Thromb Diath Haemorr</i>, 33: 26-35 Solomon G, Hilgartner M, Kutt H (1972) Coagulation defects caused by diphenylhydantoin, <i>Neurology</i>, 22: 1165-71 Ambrus JL, Schimert G, Lajos TZ, Ambrus CM, Mink IB, Lassman HB, Moore RH, Meizer J (1971) Effect of antifibrinolytic agents and estrogens on blood loss and coagulation factors during open heart surgery, <i>J Med</i>, 2: 65-81 Crowell EB Jr, Clatanoff DV, Kiekhofe W (1971) The effect of oral contraceptives on factor VIII levels, <i>Lab. Clin. Med.</i> 77: 551-7

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BIBLIOGRAFÍA
<ol style="list-style-type: none">Langdell R, Wagner R, Brinkhous K (1953) Effect of Antihemophilic Factor on One Stage Clotting Tests, <i>J. Lab. Clin. Med.</i> 41: 637 Proctor R, Rapaport S (1961) The Partial Thrombopastin Time with Kaolin, <i>Am. J. Clin. Path.</i> 36: 212 Brandt JT and Triplett DA (1981) Laboratory Monitoring of Heparin. Effect of Reagents and Instruments on the Activated Partial Thromboplastin Time, <i>AJCP</i>. 76: 530-537 Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5 Young DS <i>et al.</i> Effects of Drugs on Clinical Laboratory Tests, 3rd ed., AAC Press, Washington, D.C., 1990 Kelton JG and Hirsh J (1980) Bleeding Associated with Antithrombotic Therapy, <i>Semin Hematol</i>, 17: 259-91 Estes J and Poulin P (1974) Pharmacokinetics of Heparin, <i>Thromb Diath Haemorr</i>, 33: 26-35 Solomon G, Hilgartner M, Kutt H (1972) Coagulation defects caused by diphenylhydantoin, <i>Neurology</i>, 22: 1165-71 Ambrus JL, Schimert G, Lajos TZ, Ambrus CM, Mink IB, Lassman HB, Moore RH, Meizer J (1971) Effect of antifibrinolytic agents and estrogens on blood loss and coagulation factors during open heart surgery, <i>J Med</i>, 2: 65-81 Crowell EB Jr, Clatanoff DV, Kiekhofe W (1971) The effect of oral contraceptives on factor VIII levels, <i>Lab. Clin. Med.</i> 77: 551-7

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