REF 5502
REF 5507

Helena Biosciences Europe Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom
Tel: +44 (0)191 482 8440
Fax: +44 (0)191 482 8442
Email: info@helena-biosciences.com
Web: www.helena-biosciences.com

HL-2-1074P 2016/01 (12)

Antithrombin Xa (Chromogenic)

Instructions for use

INTENDED PURPOSE

The Antithrombin Xa (Chromogenic) kit is intended for carrying out chromogenic based haemostasis assays.

AT-III is a major inhibitor of blood coagulation, acting to inhibit plasma serine proteases including factors IXa, Xa, Xia, and thrombin. The rate of the inhibition is significantly increased in the presence of heparin. AT-III deficiency may be congenital or acquired and is associated with an increased risk for thrombosis. Deficiencies may occur in liver disease, disseminated intravascular coagulation (DIC), septicaemia, pulmonary embolism, nephrotic syndrome, stroke, and thrombophlebitis.

Most functional assays for AT-III rely on the ability of plasma to inactivate thrombin in the presence of heparin, such as the assay described by Odegaard *et al.*. More recently, it has been shown that a method based on the ability of plasma to inhibit factor Xa in the presence of heparin may be more accurate, since it eliminates interference due to heparin cofactor II[®] which does not inhibit factor Xa.

In the present two-stage method, factor Xa is added to a plasma dilution containing AT-III in the presence of excess heparin and calcium. After an initial incubation period (stage 1), residual factor Xa is determined with a factor Xa-specific chromogenic substrate (stage 2). The residual factor Xa activity is inversely proportional to the AT-III concentration. The Helena Biosciences Europe Antithrombin Xa (Chromogenic) kit is intended for the quantitative determination of antithrombin III (AT-III) activity in human citrated plasma by chromogenic assay.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

COMPOSITION

| Component | Content | Description | Preparation |
|---------------------|---|--|---|
| Factor Xa | 3 x 10 mL (REF 5502) 5 x 2 mL (REF 5507) | Contains freeze-dried bovine factor Xa. | Reconstitute each vial with: 10 mL (REF 5502) or 2 mL (REF 5507) of distilled / deionised water. Allow to stand for 10 minutes and mix well before use. |
| Factor Xa Substrate | 3 x 10 mL (REF 5502) 5 x 2 mL (REF 5507) | Contains freeze-dried CH ₃ OCO-D-CHA-Gly-Arg-pNA-AcOH. | Reconstitute each vial with: 10 mL (REF 5502) or 2 mL (REF 5507) of distilled / deionised water. Allow to stand for 10 minutes and mix well before use. |
| Sample Diluent | 4 x 10 mL (REF 5502) 5 x 3 mL (REF 5507) | Contains 5x buffer concentrate with sodium azide at 0.1% as preservative. When fully diluted, buffer contains 0.05 M Tris-HCl, 0.175 M NaCl, 7.5 mM Na ₂ EDTA and sodium heparin at pH 8.4. | Dilute 1 + 4 with distilled / deionised water and mix well. |

Each kit contains Instructions For Use.

ITEMS REQUIRED BUT NOT PROVIDED

Coagulation instrument or spectrophotometer operable at 405 nm.

Depending on Assay Method Used:

- REF 5185 Calibration Plasma
- Glacial acetic acid
- 37°C water bath or dry bath
- Laboratory timer

STORAGE, SHELF-LIFE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Factor Xa Reconstituted reagent is stable for 2 months at *2–8°C, 1 month at *17°C or 2 days at *37°C.

Factor Xa Substrate Reconstituted reagent is stable for 2 months at *2–8°C, 1 month at *17°C or 7 days at *37°C.

Sample Diluent Store diluted buffer in a tightly stoppered bottle at *2–8°C and use within one month.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma can be kept at *2–8°C or *18–24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at *37°C prior to testing. Do not keep at *37°C for more than 5 minutes.

PROCEDURE

Preparation of Standard (REF 5185), Control and Patient Dilutions:

| %AT-III | Plasma | Dilution Buffer |
|--------------------|--------------------|-----------------|
| 100% | 10 µL Standard | 990 µL |
| 50% | 500 µL of 100% Std | 500 µL |
| 25% | 500 µL of 50% Std | 500 µL |
| 12.5% | 500 µL of 25% Std | 500 µL |
| Patient or Control | 10 µL Plasma | 990 µL |

A. Semi-Micro Endpoint Method

- Add 200 µL standard / control or patient plasma dilution to a test tube.
- Incubate at *37°C for 2-4 minutes.
- Add 200 µL Factor Xa Reagent and mix.
- Incubate at *37°C for EXACTLY 1 minute.
- Add 200 µL Factor Xa Substrate and mix.
- Incubate at *37°C for EXACTLY 3 minutes.
- Add 200 µL acetic acid (50%) and mix.
- * Add 200 µL water (optional).

The yellow colour of the final reaction product is stable for at least 4 hours. Read absorbance at 405 nm in a 1 cm semi-micro cuvette against a blank prepared in the following order:

- 200 µL acetic acid.
- 200 µL standard dilution.
- 200 µL Factor Xa reagent.
- 200 µL Factor Xa Substrate.
- * 200 µL water (optional).

* Some spectrophotometers require a minimum of 1 mL volume in the cuvette.

If the patient plasma is very icteric, a second blank containing patient plasma dilution instead of standard dilution should be prepared and the absorbance subtracted from the absorbance obtained for the patient's AT-III determination. As many as ten determinations can be performed simultaneously with the same stop watch by staggering pipetting steps at five second intervals.

B. Kinetic Method

A kinetic analyser may be used to measure the initial rate of hydrolysis of the chromogenic substrate. The procedure to be used is as follows:

To a reaction cuvette:

- Add 200 µL diluted standard or patient plasma.
- Incubate at *37°C for 1 minutes.
- Add 200 µL Factor Xa reagent and mix.
- Incubate at *37°C for 1 minute.
- Add 200 µL Factor Xa Substrate.
- Measure rate of change of absorbance at 405 nm for 1 minute.

INTERPRETATION OF RESULTS

a) Assay Calibration

A commercially prepared plasma standard in which AT-III has been determined eg. Helena Biosciences Europe Calibration Plasma (REF 5185) should be used as a reference.

b) Calibration Curve

Plot the absorbance obtained with each of the AT-III calibration standards on the y-axis against % AT-III on the x-axis using linear graph paper. The line of best fit should be determined by linear regression analysis. The AT-III in plasma samples can be determined by interpolation from the calibration curve. The AT-III concentration in the patient specimen should be adjusted for the AT-III concentration in the standard:

$$\% \text{ AT-III (adjusted) } = \% \text{ AT-III (patient) } \times \% \text{ AT-III (standard) } / 100$$

LIMITATIONS

Avoid icteric, lipaemic, and haemolysed samples. To ensure reproducible results, use accurate pipetting devices and observe recommended procedures with emphasis on incubation times and temperature.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Helena Biosciences Europe supply the following controls available for use with this product:

| | |
|----------|------------------------------|
| REF 5301 | Speciality Assayed Control N |
| REF 5302 | Speciality Assayed Control A |

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range. The normal range of AT-III is 75-125% in plasma. Plasma activity levels of 30-60% can be observed in patients with hereditary AT-III deficiency^{1,2}. Several clinical conditions associated with acquired AT-III deficiency include: liver disease, DIC, nephrotic syndrome, pulmonary embolism, stroke, and thrombophlebitis. In addition, oral contraceptive use may lead to reduced AT-III levels³.

PERFORMANCE CHARACTERISTICS

Helena Biosciences Europe or their representatives have determined the following performance characteristics as a guideline. Each laboratory should establish its own performance data.

Precision

The following estimates of precision were observed using an automated analyser utilising a kinetic measurement mode:

| AT-III Level | Intra-assay precision | | Inter-assay precision | |
|--------------|-----------------------|--------|-----------------------|--------|
| | n | CV (%) | n | CV (%) |
| Normal | 20 | 1.36 | 10 | 2.34 |
| Abnormal | 19 | 4.40 | 10 | 4.99 |

Sensitivity

This assay is designed to give a linear standard curve for AT-III levels between 7.5-150%. Samples over 150% should be diluted in Dilution Buffer and re-assayed.

Specificity

The specificity of the assay system has been established in studies employing plasma that has been selectively depleted of AT-III followed by the addition of purified AT-III to achieve various AT-III concentrations.

Accuracy

Comparison of the Helena Biosciences Europe Antithrombin Xa (Chromogenic) kit with another commercially available kit for AT-III determination on 62 samples ranging from 30-130% AT-III activity, the following results were obtained:

$$y = 0.9665x + 0.0629 \quad r^2 = 0.8893$$

BIBLIOGRAPHY

- Hirsh J *et al* (1989) Congenital antithrombin III deficiency. Incidence and clinical features, *Am J Med*, **87**(suppl. 3B): 34S-38S
- Panicucci F *et al* (1981) Antithrombin III, heparin cofactor and antifactor Xa in relation to age, sex and pathological condition, *Haemostasis*, **9**: 297-302
- Conlan MG *et al* (1994) Antithrombin: Antithrombin III: Associations with age, race, sex and cardiovascular disease risk factors, *Thromb. Haemost.*, **72**(4): 551-556
- Odegaard OR *et al* (1975) Heparin cofactor activity measured with an amidolytic method, *Thromb. Res.*, **6**: 287-294
- Demers C *et al* (1993) An antithrombin III assay based on factor Xa inhibition provides a more reliable test to identify congenital antithrombin III deficiency than an assay based on thrombin inhibition, *Thromb. Haemost.*, **69**: 231-235
- Tolefsen DM, Blank MK (1981) Detection of a new heparin-dependent inhibitor of thrombin in human plasma, *J. Clin. Invest.*, **68**: 589-596
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5

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