
Auto Red D-Dimer 700

REF 5501

	Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom
	Tel: +44 (0)191 482 8440
	Fax: +44 (0)191 482 8442
	Email: info@helena-biosciences.com
	Web: www.helena-biosciences.com

HL-2-1790P 2015/10 (5)

Auto Red D-Dimer 700 Instructions for use

INTENDED PURPOSE

The Auto Red D-Dimer 700 kit is intended for carrying out immunoturbidimetric based haemostasis assays.

D-dimer containing moieties are formed by plasmin degradation of factor XIIIa cross-linked fibrin. Elevated levels of D-dimer are found in clinical conditions such as deep vein thrombosis (DVT), pulmonary embolism (PE) and disseminated intravascular coagulation (DIC)¹⁻³. Laboratory measurements of fibrin degradation products, including D-dimer, have significance in the initial assessment of these conditions.

Auto Red D-Dimer 700 is a turbidimetric assay that utilises antibody coated latex particles. In the presence of D-dimer, the particles aggregate and turbidity increases. The increase in scattered light is proportional to the amount of D-dimer in the sample. The latex particles are coated with a monoclonal antibody that reacts with fibrin D-dimer or fragment D of fibrin. The antibody has no cross reactivity with fibrinogen⁴. This allows for the determination of D-dimer in human plasma. Auto Red D-Dimer 700 is an immunoturbidimetric assay used for the quantitative determination of the fibrin degradation products that contain D-dimer in human plasma.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of:
Hepatitis B Antigen (HbsAg)
HIV 1 antibody
HIV 2 antibody
HCV antibody
However they should be handled with the same precautions as a human patient sample.

COMPOSITION

Important: The reagents are lot-specific. Lots are not interchangeable.

Component	Content	Description	Preparation
D-Dimer Red Latex	4 x 4 mL (REF 5501)	Latex particles coated with anti-D-dimer monoclonal antibody.	Ready to use. The latex may sediment during storage. Mix thoroughly before use.
D-Dimer Red Buffer	4 x 7 mL (REF 5501)	Containing buffer and preservatives.	Ready to use.
D-Dimer Diluent	2 x 7 mL (REF 5501)	Containing buffer and preservatives.	Ready to use.
D-Dimer Calibrator	2 x 1 mL (REF 5501)	Lyophilised human plasma enriched with D-dimer.	Reconstitute each vial of the appropriate calibrator with 1 mL of D-Dimer Diluent. Swirl gently. Allow to stand for 15 minutes for complete dissolution and mix well before use.

Each kit contains instructions for use.

Each kit contains lot specific reference values insert.

ITEMS REQUIRED BUT NOT PROVIDED

Analysér operable at 600-900 nm.
Helena Biosciences Europe supplies the following controls available for use with this product:

REF 5509 D-Dimer Control H/L

The reagents contained in these kits are also available as separate items:

REF 32097SA D-Dimer Calibrator
REF 32096SA D-Dimer Diluent

STORAGE, SHELF-LIFE AND STABILITY

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

D-Dimer Red Latex Once opened, the reagent is stable for 4 weeks at *2 –*8°C or 2 weeks at *20°C.

D-Dimer Red Buffer Once opened, the reagent is stable for 4 weeks at *2 –*8°C or 2 weeks at *20°C.

D-Dimer Diluent Store at *2 –*8°C and use within 4 weeks of opening.

D-Dimer Calibrator Once reconstituted, the reagent is stable for 12 hours at *4 –*25°C.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes⁵. Plasma should be kept at *2 –*8°C or *18 –*24°C. Testing should be completed within 24 hours of sample collection, or plasma can be stored frozen at -20°C or -70°C for 24 months⁶. Thaw quickly at *37°C prior to testing. Do not keep at *37°C for more than 5 minutes.

PROCEDURE

Preparation of Standard Curve

Use the lot-specific value for D-Dimer Calibrator to determine the exact D-dimer concentration in each standard dilution. Users must construct a standard curve each time a new kit lot is used, after any instrument change or service, or every 6 months, and if D-Dimer Control H/L (REF 5509) is assayed out of range.

Automated Method

Refer to the appropriate instrument operator manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application guides.

INTERPRETATION OF RESULTS

The test should be used in conjunction with clinical observations and results of other laboratory tests. Elevated levels are found in patients with confirmed DVT, PE, DIC, and trauma¹⁻³. D-dimer levels rise during pregnancy and high levels are associated with complications⁷. The results may be reported either in D-dimer units or in fibrinogen equivalent units (FEU); 1 ng/mL of D-dimer is equivalent to approximately 2 ng/mL of FEU, however, a more accurate conversion factor, originating from the Fibrinogen/D-dimer weight ratio of 340 kDa/195 kDa, would be 1.74⁸.

LIMITATIONS

Presence of rheumatoid arthritis factor may result in false-positive results (influence not quantified). Results from patients with heterophilic antibody should be interpreted with caution since this test kit contains mouse antibodies and interference may occur resulting in falsely elevated or decreased values. Turbid or opalescent plasma may cause erratic results and should be interpreted with caution.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. It is recommended that the D-dimer low control and D-dimer high control are assayed at regular intervals in order to ensure consistent assay results. If the D-Dimer Control H/L result deviates from the D-dimer concentration given in the lot-specific Instruction for Use, a new standard curve should be constructed.

REFERENCE VALUES

The concentration of D-dimer in any given specimen may differ from the concentration determined using D-dimer assays from different manufacturers. Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives using a photo-optical coagulation instrument. The user should establish product performance characteristics for the specific instrumentation used.

Reproducibility	Intra-assay precision	Inter-assay precision			
<i>Sample</i>	<i>n</i>	<i>D-dimer (ng/mL)</i>	<i>CV (%)</i>	<i>D-dimer (ng/mL)</i>	<i>CV (%)</i>
Medium	42	1556	2,1	1541	1,1
Low	42	304	5,0	308	1,8

Auto Red D-Dimer 700 is insensitive to the following substances: haemoglobin (up to 10 g/L), bilirubin (up to 0.5 g/L), triglycerides (up to 20 g/L), low molecular weight heparin (up to 100 U/mL), and non-fractionated heparin (up to 100 U/mL).

BIBLIOGRAPHY

- Declerck P *et al.* (1987) Fibrinolytic response and fibrin fragment D-dimer levels in patients with deep vein thrombosis, *Thrombosis and Haemostasis*, **58**:1024-1029
- Lindahl T *et al.* (1998) Clinical evaluation of a diagnostic strategy for deep venous thrombosis with exclusion by low plasma levels of fibrin degradation product D-dimer, *Scand. J. Clin. Lab. Invest*, **58**:307-316
- Hansson PO *et al.* (1994) Can laboratory testing improve screening strategies for deep vein thrombosis at an emergency unit? *J. Intern. Med*, **235**:143-151
- Holvoet P *et al.* (1989) Binding properties of monoclonal antibodies against human fragment D-dimer of cross-linked fibrin to human plasma clots in an in vivo model in rabbits, *Thrombosis and Haemostasis*, **61**:307-313
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Clinical and Laboratory Standards Institute (2011) Quantitative D-dimer for the Exclusion of Venous Thromboembolic Disease; Approved Guideline, CLSI: H59-A
- Ballegeer V *et al.* (1987) Fibrinolytic response to venous occlusion and fibrin fragment D-dimer levels in normal and complicated pregnancy, *Thrombosis and Haemostasis*, **58**:1030-1032
- Edlund B, Nilsson TK (2006) A proposed stoichiometrical calibration procedure to achieve transferability of D-dimer measurements and to characterize the performance of different methods, *Clin Biochem*, **39**(2):137-142

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