Clauss Fibrinogen 100

REF 5376
REF 5376H

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom
Tel: +44 (0)191 482 8440
Fax: +44 (0)191 482 8442
Email: info@helena-biosciences.com
Web: www.helena-biosciences.com

HL-2-0592P 2016/01 (11)

Clauss Fibrinogen 100 Instructions for use

INTENDED PURPOSE

The Clauss Fibrinogen 100 kit is intended for carrying out clot based haemostasis assays.

Clauss¹ developed a simple method for the quantitative determination of fibrinogen by measuring the clotting time of dilute plasma after the addition of Thrombin (>30 NIH units/mL). This clot time is proportional to the fibrinogen concentration. Levels of fibrinogen can increase as a result of inflammation, pregnancy or oral contraceptive use². Decreased levels can be found in certain states such as liver disease and DIC. Congenital deficiencies include afibrinogenaemia (no detectable fibrinogen), hypofibrinogenaemia (<1 mg/mL) and dysfibrinogenaemia (abnormal fibrinogen molecule). Clauss Fibrinogen 100 is intended for the quantitative determination of fibrinogen in human plasma.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of:

- Hepatitis B Antigen (HbsAg)
- HIV 1 antibody HIV 2 antibody
- HCV antibody

However they should be handled with the same precautions as a human patient sample.

COMPOSITION

Component	Content	Description	Preparation
Thrombin: 100 NIH/mL	5 x 2 mL (REF 5376) 5 x 1 mL (REF 5376H)	Contains approximately 200 (REF 5376) / 100 (REF 5376H) units bovine thrombin with stabilisers.	Reconstitute each vial with: 2 mL (REF 5376) 1 mL (REF 5376H) of purified water or kaolin suspension depending on manual or instrument specific methodology. Take care when pipetting to avoid contamination. Swirl gently and allow to stand for 15 minutes. Mix well immediately before use. Do not shake.
Fibrinogen Calibrator	2 x 1 mL (REF 5376) 1 x 1 mL (REF 5376H)	Contains 1 mL of lyophilised normal human plasma assayed for fibrinogen levels.	Reconstitute each vial with exactly 1 mL of purified water. Swirl gently and allow to stand for 10 minutes. Mix gently before use. Do not shake.
Owren's Buffer	2 x 25 mL (REF 5376) 1 x 25 mL (REF 5376H)	Contains 25 mL of buffer which contains barbital and sodium chloride and sodium azide as preservative.	The buffer is ready for use as packaged.
Kaolin Suspension: 0.5 g/L	2 x 5 mL (REF 5376) 1 x 5 mL (REF 5376H)	Contains 5 mL of kaolin suspension at 0.5 g/L.	Mix vigorously immediately before use. Ensure no sedimentation is visible.

Each kit contains instructions for use.

Each kit contains lot specific reference values insert.

ITEMS REQUIRED BUT NOT PROVIDED

Helena Biosciences Europe fibrinogen assay reagents are suitable for use on any coagulation analyser capable of performing Clauss Fibrinogen assays (not for use with the AC-4, C-1, C-2, C-4). Refer to the Instrument Operators' Manual and/or application note for appropriate instructions.

The reagents contained in this kit are also available as separate items:

REF 5374	Clauss Fibrinogen (Thrombin only)
REF 5378	Clauss Fibrinogen (Thrombin only)
REF 5379	Fibrinogen Calibrator
REF 5375	Owren's Buffer
REF 53765	Kaolin Suspension

STORAGE, SHELF-LIFE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Thrombin: 100 NIH/mL	Once reconstituted, the reagent is stable for 8 hours at *15 –*30°C, 1 week at *2 –*8°C or 1 month at -20°C.
Fibrinogen Calibrator	Once reconstituted, the reagent is stable for 4 hours at *2 –*8°C.
Owren's Buffer	Store at *2 –*8°C once opened.
Kaolin Suspension	Store at *2 –*8°C once opened.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at *2 –*8°C or *18 –*24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at *37°C prior to testing. Do not keep at *37°C for more than 5 minutes³.

PROCEDURE

Manual Method

Prepare all reagents as instructed under "composition".

- Standard Curve Preparation:
 - A new standard curve should be used when reagent lot numbers change or if expected values in quality control testing change. Prepare the following dilutions in Owren's Buffer (mix without shaking):

Tube	Dilution	Fibrinogen Calibrator (mL)	Buffer (mL)
1	1 + 4	0.2	0.8
2	1 + 9	0.1	0.9
3	1 + 19	0.1	1.9
4	1 + 29	0.1	2.9
5	1 + 39	0.1	3.9

- Patient Sample Preparation:
 - Prepare 1+9 dilutions of the patient plasma or control plasma in Owren's Buffer.
 - Mix without shaking.
- Testing:
 - Perform all tests in duplicate:
 - Pipette 0.2 mL of standard, patient or control dilution into a reaction tube and incubate at *37°C for 2 minutes.
 - Add 0.1 mL of thrombin (*15 –*30°C).
 - Determine the clot time to the nearest 0.1 seconds.
 - Plot mean Standard Clot Time (Y-axis) versus Fibrinogen Level (X-axis) on the Fibrinogen graph paper. A straight line should be obtained. Assign the Fibrinogen Calibrator reference value to the 1+9 dilution to allow direct interpolation of patient and control values from the standard curve.

Automated Method

Refer to the appropriate instrument Operator Manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application guides.

INTERPRETATION OF RESULTS

Expected values for fibrinogen in healthy adults are 150-350 mg/dL (1.5-3.5 g/L)^{4,5}.

LIMITATIONS

Heparin levels >0.6 units/mL and fibrinolytic degradation products >100 mg/mL may cause falsely low fibrinogen quantitation. If values fall outside the standard curve values for the patient samples, re-assay using an appropriate dilution to bring values into the standard range.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. Helena Biosciences Europe supply the following controls available for use with this product:

REF 5301	Speciality Assayed Control N
REF 5302	Speciality Assayed Control A
REF 5186	Routine Control N
REF 5187	Routine Control A
REF 5183	Routine Control SA

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range.

PERFORMANCE CHARACTERISTICS

Helena Biosciences Europe or their representatives have determined the following performance characteristics as a guideline. Each laboratory should establish its own performance data. The manual Clauss Fibrinogen is designed to give a linear calibration from 1.5 - 6.5 g/L.

Reproducibility			Inter-assay precision		
Intra-assay precision			Inter-assay precision		
<i>Fibrinogen (g/L)</i>	<i>n</i>	<i>CV (%)</i>	<i>Fibrinogen (g/L)</i>	<i>n</i>	<i>CV (%)</i>
1.24	5	4.88	1.02	10	6.17
3.01	5	3.58	3.62	10	3.75

BIBLIOGRAPHY

- Clauss A (1957) Gerinnungsphysiologische Schnell-methode zur bestimmung des Fibrinogens, *Acta Haematol*, **17**:237-246
- Shaw TS (1977) Assays for Fibrinogen and its Derivatives CRC, *Crit. Rev. Clin. Lab. Sci.* **8**:145-192.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Scully RE et al. (1980) Normal Reference Laboratory Values, *New England Journal Medicine*, **302**(37):37-48.
- Okuno T, Selenko V (1972) Plasma fibrinogen determination by automated thrombin time, *American Journal of Medical Technology*, **38**(6):196-201

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BIBLIOGRAFIA

- Clauss A (1957) Gerinnungsphysiologische Schnell-methode zur bestimmung des Fibrinogens, *Acta Haematol*, **17**:237-246
- Shaw TS (1977) Assays for Fibrinogen and its Derivatives CRC, *Crit. Rev. Clin. Lab. Sci.* **8**:145-192.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Scully RE *et al.* (1980) Normal Reference Laboratory Values, *New England Journal Medicine*, **302**(37):37-48.
- Okuno T, Selenko V (1972) Plasma fibrinogen determination by automated thrombin time, *American Journal of Medical Technology*, **38**(6):196-201

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