Clauss Fibrinogen 50

REF 5556R
REF 5556HR

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate,
Gateshead, Tyne and Wear, NE11 0SD, United Kingdom

Tel: +44 (0)191 482 8440

Fax: +44 (0)191 482 8442

Email: info@helena-biosciences.com

Web: www.helena-biosciences.com

HL-2-2739P 2013/07 (1)

Clauss Fibrinogen 50 Instructions for use

en

INTENDED USE

For the quantitative determination of fibrinogen, based on the Clauss method, in human citrated plasma on Helena Biosciences Europe AC-4, C-1, C-2, C-4 Coagulation Systems.

Clauss Fibrinogen 50 is intended for the quantitative determination of fibrinogen in human plasma. Clauss¹ developed a simple method for the quantitative determination of fibrinogen by measuring the clotting time of dilute plasma after the addition of Thrombin (>30 NIH units/mL). This clot time is proportional to the fibrinogen concentration. Levels of fibrinogen can increase as a result of inflammation, pregnancy or oral contraceptive use². Decreased levels can be found in certain states such as liver disease and DIC. Congenital deficiencies include afibrinogenemia (no detectable fibrinogen), hypofibrinogenemia (<1 mg/mL) and dysfibrinogenemia (abnormal fibrinogen molecule).

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only - DO NOT INGEST. Wear gloves when handling all kit components. Refer to the product safety data sheets for risk and safety phrases and disposal information. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 and 2 antibody and HCV antibody; however they should be handled with the same precautions as a human plasma sample.

COMPOSITION

Component	Content	Description	Preparation
Thrombin: 50 NIH/mL	5 x 4 mL (REF 5556R) 5 x 2 mL (REF 5556HR)	Contains approximately 200 (REF 5556R) / 100 (REF 5556HR) units bovine thrombin with stabilisers.	Reconstitute each vial with: 4 mL (REF 5556R) 2 mL (REF 5556HR) of purified water. Take care when pipetting to avoid contamination. Swirl gently and allow to stand for 15 minutes. Mix well immediately before use. Do not shake.
Fibrinogen Calibrator	2 x 1 mL (REF 5556R) 1 x 1 mL (REF 5556HR)	Contains 1 mL of lyophilised normal human plasma assayed for fibrinogen levels.	Reconstitute each vial with exactly 1 mL of purified water. Swirl gently and allow to stand for 10 minutes. Mix gently before use. Do not shake.
Imidazole Buffer	2 x 25 mL (REF 5556R) 1 x 25 mL (REF 5556HR)	Each bottle contains an aqueous solution of 0.05M Imidazole, sodium chloride, and sodium azide (<0.1%) as a preservative pH 7.3 ±0.10.	The buffer is ready for use as packaged.

Each kit contains instructions for use.

Each kit contains lot specific reference values insert.

ITEMS REQUIRED BUT NOT PROVIDED

Refer to the Instrument Operator Manual and/or application note for appropriate instructions. The reagents contained in this kit are also available as separate items:

REF 5374	Clauss Fibrinogen (Thrombin only)
REF 5378	Clauss Fibrinogen (Thrombin only)
REF 5379	Fibrinogen Calibrator
REF 5375R	Imidazole Buffer

STORAGE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Thrombin: 50 NIH/mL Once reconstituted, the reagent is stable for 8 hours at +15 –+30°C, 1 week at +2 –+8°C or 1 month at -20°C.

Fibrinogen Calibrator Once reconstituted, the reagent is stable for 4 hours at +2 –+8°C.

Imidazole Buffer Store at +2 –+8°C once opened.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at +2 –+8°C or +18 –+24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at +37°C prior to testing. Do not keep at +37°C for more than 5 minutes³.

PROCEDURE

Refer to the appropriate instrument Operator Manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application guides.

INTERPRETATION OF RESULTS

Expected values for fibrinogen in healthy adults are 150-350 mg/dL (1.5-3.5 g/L)^{4,5}.

LIMITATIONS

Heparin levels >0.6 units/mL and fibrinolytic degradation products >100 mg/mL may cause falsely low fibrinogen quantitation. If values fall outside the standard curve values for the patient samples, re-assay using an appropriate dilution to bring values into the standard range.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. Helena Biosciences Europe supply the following controls available for use with this product:

REF 5301	Speciality Assayed Control N
REF 5302	Speciality Assayed Control A
REF 5186	Routine Control N
REF 5187	Routine Control A
REF 5183	Routine Control SA

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range.

PERFORMANCE CHARACTERISTICS

Helena Biosciences Europe or their representatives have determined the following performance characteristics as a guideline. Each laboratory should establish its own performance data. The Clauss Fibrinogen assay is designed to give a linear calibration from 1.5 - 6.5 g/L.

Reproducibility

Intra-assay precision

Fibrinogen (g/L)	n	CV (%)
1.2	21	5.7
3.2	21	5.9
6.3	21	1.8

REFERENCES

- Clauss A (1957) Gerinnungsphysiologische Schnell-methode zur bestimmung des Fibrinogens, *Acta Haematol*, **17**:237-246
- Shaw TS (1977) Assays for Fibrinogen and its Derivatives CRC, *Crit. Rev. Clin. Lab. Sci*, **8**:145-192.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Scully RE *et al.* (1980) Normal Reference Laboratory Values, *New England Journal Medicine*, **302**(37):37-48.
- Okuno T, Selenko V (1972) Plasma fibrinogen determination by automated thrombin time, *American Journal of Medical Technology*, **38**(6):196-201

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