D-Dimer Control

helenabiosciences Europe

REF 5509 D-Dimer Control H/L

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 0SD, United Kingdom

Tel: +44 (0)191 482 8440

Fax: +44 (0)191 482 8442

Email: info@helenabiosciences.com

Web: www.helenabiosciences.com

HL-2-1188P 2015/10 (10)

D-Dimer Control Instructions for use

en

INTENDED PURPOSE

The D-Dimer Control H/L kit is intended for use as a quality control material.

D-dimer containing moieties are formed by plasmin degradation of Factor XIIIa cross-linked fibrin. Elevated levels of D-dimer are found in clinical conditions such as deep vein thrombosis (DVT), pulmonary embolism (PE) and disseminated intravascular coagulation (DIC). D-dimer levels rise during pregnancy and high levels are associated with complications. The negative predictive value of D-dimer for thrombosis is high, although a negative D-dimer does not completely rule out thrombosis. The Helena Biosciences Europe D-Dimer Control H/L kit contains plasmas with low and high levels of D-dimer. These plasmas are for use with the Helena Biosciences Europe latex-enhanced immunoassay methods for D-dimer and should be tested before patient samples to ensure the satisfactory performance of the assay system.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of: Hepatitis B Antigen (HbsAg), HIV 1 antibody, HIV 2 antibody.

However they should be handled with the same precautions as a human patient sample.

COMPOSITION

Component	Content	Description	Preparation
D-Dimer Control - L	5 x 1 mL (REF 5509)	Each vial contains lyophilised D-dimer enriched human plasma.	Reconstitute each vial with 1 mL of saline solution. Allow to stand for 5 minutes and mix well before use. Do not shake.
D-Dimer Control - H	5 x 1 mL (REF 5509)	Each vial contains lyophilised D-dimer enriched human plasma.	Reconstitute each vial with 1 mL of saline solution. Allow to stand for 5 minutes and mix well before use. Do not shake.

Each kit contains instructions for use.

Each kit contains lot specific reference values insert.

ITEMS REQUIRED BUT NOT PROVIDED

REF 5551 / 5552 / 5554	Auto Blue D-Dimer 400
REF 5501 / 5601	Auto Red D-Dimer 700

STORAGE, SHELF-LIFE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label. Reconstituted vials are stable for 5 days at 2 – 8 °C, or 3 months at -20°C. Do not freeze / thaw more than once.

SAMPLE COLLECTION AND PREPARATION

Not applicable.

PROCEDURE

Each control should be treated in the same manner as the unknown specimen in accordance with the instructions outlined in each particular test protocol (REF 5551, 5552, 5554, 5501, 5601).

INTERPRETATION OF RESULTS

Lot specific expected values are provided with each pack of controls. Recalibration may be required if results fall outside quoted ranges.

LIMITATIONS

The values listed are for use with Helena Biosciences Europe D-dimer procedures only. The use of other D-dimer procedures and products may lead to erroneous results.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives using an opto-mechanical coagulation instrument. Each laboratory should establish its own performance data.

Reproducibility

Sample	n	Intra-assay precision CV (%)	Inter-assay precision CV (%)
D-Dimer Control - L	35	3.3	2.5
D-Dimer Control - H	35	6.4	0.9

BIBLIOGRAPHY

1. Elms MJ *et al* (1986) Rapid Detection of Cross-Linked Fibrin Degradation Products in Plasma Using Monoclonal Antibody-Coated Latex Particles, *J. Clin. Pathol.* **35**: 360-64.
2. Declercq PV *et al* (1987) Fibrinolytic Response and Fibrin Fragment D-dimer in Patients With Deep Vein Thrombosis, *Thromb. Haemost.* **58**: 1025-9.
3. Balleger V *et al* (1987) Fibrinolytic Response to Venous Occlusion and Fibrin Fragment D-dimer Levels in Normal and Complicated Pregnancy, *Thromb. Haemost.* **58**: 1030-2.
4. Hansson PO *et al* (1994) Can Laboratory Testing Improve Screening Strategies for Deep Vein Thrombosis at an Emergency Unit?, *J. Intern. Med.* **235**: 143-51.

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