REF 5484

| | |
|--|--|
| CE IVD | |
|  | Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom |
| | Tel: +44 (0)191 482 8440 |
| | Fax: +44 (0)191 482 8442 |
| | Email: info@helena-biosciences.com |
| | Web: www.helena-biosciences.com |

HL-2-0723P 2015/10 (7)

| | |
|-------------------------|-----------|
| DRVVT Screen | en |
| Instructions for use | |
| INTENDED PURPOSE | |

The DRVVT Screen kit is intended for carrying out clot based haemostasis assays.

Lupus Anticoagulants (LA) are antibodies of the IgG and IgM type which are directed against a variety of anionic phospholipids. The presence of LA in plasma is increasingly associated with a variety of haemostatic problems such as recurrent fetal loss, thrombocytopenia, unexplained thrombosis and neurological disorders. LA prolongs phospholipid dependant *in vitro* clotting assays such as the activated partial thromboplastin time (aPTT). The Helena Biosciences Europe DRVVT Screen kit is intended for the qualitative determination of LA in human plasma. Russell's Viper venom directly activates Factor X to Factor Xa in the presence of phospholipid and calcium, leading to detectable clot formation in plasma. The DRVVT Screen kit is more sensitive for LA than the aPTT. The DRVVT Screen kit is intended to be used in conjunction with the DRVVT Confirm kit (REF 5485).

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

| Component | Content | Description | Preparation |
|--------------|-----------|--|--|
| DRVVT Screen | 10 x 2 mL | Each vial contains a proprietary mixture of Russell's Viper venom co-lyophilised with calcium chloride and phospholipid. | Reconstitute each vial with 2 mL of distilled or deionised water. Allow to stand for 10 minutes and mix well before use. Do not shake. |

Each kit contains Instructions For Use.

ITEMS REQUIRED BUT NOT PROVIDED

OPTIONAL ITEMS:

| | |
|----------|-----------------------|
| REF 5485 | DRVVT Confirm |
| REF 5486 | LA Positive Control S |
| REF 5186 | Routine Control N |

STORAGE, SHELF-LIFE AND STABILITY

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

DRVVT Reconstituted vials are stable for 24 hours at *15 –*30°C, 5 days at *2 –*8°C or 2 weeks at Screen -20°C. The reagent should be frozen in plastic test tubes and thawed at *37°C before use.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at *2 –*8°C or *18 –*24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at *37°C prior to testing. Do not keep at *37°C for more than 5 minutes¹. If freezing, double centrifugation of the sample is recommended to ensure that the sample is platelet poor. Transfer the plasma following the initial centrifugation to a non-activating plastic tube using a plastic pipette, then re-centrifuge the plasma for an additional 10 minutes at a higher speed (>2500 x g). When aliquoting to a secondary tube, take care to not include the residual platelets that may have collected at the bottom of the centrifuge tube².

PROCEDURE

Manual Method

- Pre-warm sufficient reconstituted reagent to *37°C.
- Pipette 0.2 mL of patient or control plasma into a reaction tube. Incubate at *37°C for 2 minutes.
- Add 0.2 mL of pre-warmed DRVVT Screen Reagent and start a timer.
- Measure the clot formation time to the nearest 0.1 seconds.
- Calculate the normalised 'DRVVT Screen' ratio as:

Patient DRVVT Screen Clot Time / Mean Normal DRVVT Screen Clot Time

Automated Method

Refer to the appropriate instrument operator manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application guides.

INTERPRETATION OF RESULTS

Results are best expressed as a normalised ratio relative to the mean normal clot time obtained by each laboratory. It is recommended that like for like sample types are used when calculating a normalised ratio. Both DRVVT Screen and DRVVT Confirm results can be 'normalised' in this way, reducing the effect of instrument variability and potentially improving discrimination between weak positive LA and normal samples. Results of the mixing tests can be treated in the same way.

If the clot time of the patient sample is greater than 3 standard deviations above the mean of the normal range, a lupus anticoagulant may be present. In this case, the plasma should be re-tested after mixing 1:1 with Routine Control N (REF 5186) as well as testing with the DRVVT Confirm kit (REF 5485). If the DRVVT Screen clotting time of the patient plasma mixed 1:1 with Routine Control N is still greater than 3 standard deviations above the mean of the normal range, a lupus anticoagulant may be present. If the DRVVT Screen clotting time of the patient plasma mixed 1:1 with Routine Control N is corrected to within the normal range, a factor deficiency (II, V or X) is most likely.

The Scientific and Standardisation Sub-Committee for the Standardisation of Lupus Anticoagulants of the International Society of Thrombosis and Haemostasis has recommended that the diagnosis of lupus anticoagulant be made when the DRVVT of a test plasma mixed with normal plasma is greater than 3 standard deviations from the mean normal (non-LA) plasma DRVVT time.

The use of the DRVVT Confirm kit allows discrimination between LA, factor deficiency and other inhibitors.

LIMITATIONS

Plasma deficiencies of Factors II, V or X may lead to abnormal results in neat plasma. Mixing studies should correct this.

Plasma from patients with the following may give abnormal results when the plasma is tested neat, and these samples may not correct in mixing studies: heparin (>1 U/mL), oral anticoagulants, disseminated intravascular coagulation (DIC).

Care must be taken to remove residual platelets from plasma by filtration or centrifugation, as platelet derived phospholipid can interfere with the test.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. Helena Biosciences Europe supplies the following controls available for use with this product:

| | |
|----------|-----------------------|
| REF 5486 | LA Positive Control S |
| REF 5186 | Routine Control N |

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges. The normal reference range (mean ± 3SDs) determined at Helena Biosciences Europe for the DRVVT Screen test was 33.7 ± 8.1 seconds (range 25.6 - 41.8 seconds).

PERFORMANCE CHARACTERISTICS

Each laboratory should establish its own performance data. Within run and between run CVs are expected to be <5%.

BIBLIOGRAPHY

- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Pengo *V et al* (2009) Update of the guidelines for lupus anticoagulant detection, *J Thromb Haemost*, **7**: 1737-40

| | |
|---------------------|-----------|
| DRVVT Screen | fr |
| Fiche technique | |

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