Epinephrine

helena
Biosciences Europe

REF 5367

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 0SD, United Kingdom

Tel: +44 (0)191 482 8440

Fax: +44 (0)191 482 8442

Email: info@helena-biosciences.com

Web: www.helena-biosciences.com

HL-2-0500P 2015/10 (13)

Epinephrine Instructions for use

en

INTENDED PURPOSE

The Epinephrine kit is intended for use in platelet aggregation studies.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate advice and precautionary statements where applicable. Dispose of components in accordance with local regulations.

COMPOSITION

Composition	Content	Description	Preparation
Epinephrine: 3 mM	2 x 1 mL	Each vial contains 1 mg of epinephrine bitartrate, which is equivalent to 548 µg (3 µmoles) epinephrine. The reagent should be a glassy solid or uniform white powder before reconstitution.	Reconstitute each vial with 1 mL of distilled or deionised water. Mix gently before use ensuring complete dissolution. Following reconstitution each vial contains a clear colourless solution of epinephrine at 3 mM.

Each kit contains Instructions For Use.

ITEMS REQUIRED BUT NOT PROVIDED

Platelet Aggregometer
0.9% Saline

STORAGE, SHELF-LIFE AND STABILITY

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Epinephrine: 3 mM The reconstituted reagent is stable for 1 week when stored at *2 –8°C or 8 weeks at -20°C when flash frozen.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate platelet-rich plasma after centrifugation at 170 x g for 10 minutes (*18 –25°C). Separate platelet-poor plasma after centrifugation at 1500 x g for 15 minutes. Dilute the platelet-rich plasma with platelet-poor plasma to obtain a platelet count of 200-250 x 10⁹/L. Plasma should be kept at *18 –25°C. Testing should be completed within 4 hours of sample collection.

PROCEDURE

No dilution is necessary prior to use, although if dilution is required it should be performed in saline and used immediately. Refer to the appropriate instrument operator manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application notes.

INTERPRETATION OF RESULTS

Several disease states may be identified by interpretation of patients' plasma response to each agonist. Aggregation response should be determined by evaluating maximum aggregation, final aggregation, lag phase, primary slope, secondary slope and time to maximum aggregation values for each agonist tested. Please refer to laboratory standard operating procedure for comprehensive result interpretation in relation to each disease state. Each laboratory should determine its own reference intervals for each agonist¹.

LIMITATIONS

Avoid icteric, lipaemic, and haemolysed samples.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal controls should be tested prior to each batch of patient samples to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

PERFORMANCE CHARACTERISTICS

Within run CVs for Helena Biosciences Europe aggregation reagents using the PACKS-4 were:

Reproducibility

Intra-assay precision

Sample	n	Maximum Aggregation	CV (%)
Normal Pooled Platelet-Rich Plasma	10	91.37	2.08

BIBLIOGRAPHY

1. CLSI (2008) Platelet function testing by aggregometry; Approved guideline. H58-A Vol 28 No.31.

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