Factor Deficient Plasmas

(APTT Based - Factors VIII, IX, XI, XII)

- REF 5193 Factor VIII Deficient Plasma (Congenital)
- REF 5194 Factor IX Deficient Plasma (Congenital)
- REF 5196 Factor XI Deficient Plasma (Congenital)
- REF 5197 Factor XII Deficient Plasma (Congenital)
- REF 5793 Factor VIII Deficient Plasma (Immunodepleted)
- REF 5794 Factor IX Deficient Plasma (Immunodepleted)
- REF 5796 Factor XI Deficient Plasma (Immunodepleted)
- REF 5797 Factor XII Deficient Plasma (Immunodepleted)

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom
Tel: +44 (0)191 482 8440
Fax: +44 (0)191 482 8442
Email: info@helena-biosciences.com
Web: www.helena-biosciences.com

HL-2-0600P 2016/01 (12)

Factor Deficient Plasmas

(APTT Based - Factors VIII, IX, XI, XII)

Instructions for use

INTENDED PURPOSE

The Factor Deficient Plasmas are intended for carrying out clot based haemostasis assays.

The Factor Deficient Plasmas are intended for the quantitative determination of the respective factor in patients suspected of having a congenital or acquired deficiency of this coagulation protein. Quantitative measurement of individual coagulation factors by the one-stage method requires substrate plasma lacking the factor to be measured. A dilution of the test plasma is mixed with the factor deficient plasma and the clot time of the mixture determined. The degree of clot time correction with the patient plasma is compared to the correction with a reference material, allowing the % activity of the patient plasma to be determined¹. The Helena Biosciences Europe Factor Deficient Plasmas can be used on any instrument capable of performing APTT-based factor assay testing. Refer to the instrument Operators' Manual for appropriate instructions.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of:
Hepatitis B Antigen (HbsAg)
HIV 1 antibody
HIV 2 antibody
HCV antibody
However they should be handled with the same precautions as a human patient sample.

COMPOSITION

Component	Content	Description	Preparation
Factor Deficient Plasmas	10 x 1 mL	All the Factor Deficient Plasmas listed above are derived from human plasma and contain less than 1% residual factor activity.	Reconstitute each vial with 1 mL of purified water. Swirl gently and allow to stand for 15 minutes. Mix well before use (do not shake).

Each kit contains Instructions For Use.

ITEMS REQUIRED BUT NOT PROVIDED

REF 5185	Calibration Plasma
REF 5375	Owren's Buffer
REF 5562 / 5558 / 5559 / 5562SLQ / 5558SLQ / 5559SLQ	APTT Si L Minus

STORAGE, SHELF-LIFE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label. Once reconstituted, the reagent is stable for 8 hours when kept at *2 –*8°C. The lyophilised product should appear as a dry, straw coloured plug or pieces. Any deviation from this appearance may indicate signs of product deterioration.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at *2 –*8°C or *18 –*24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at *37°C prior to testing. Do not keep at *37°C for more than 5 minutes*.

PROCEDURE

Manual Method

Prepare all reagents as instructed with each pack. Pre-warm both the APTT Si L Minus (REF 5562 / 5558 / 5559/ 5562SLQ / 5558SLQ / 5559SLQ) reagent and calcium chloride solution to *37°C.

- Standard Curve Preparation:
 - Prepare the following dilutions in Owren's Buffer:

Tube	Calibration Plasma (mL) (REF 5185)	Owren's Buffer (mL) (REF 5375)	Activity (%)
1	0,1	0,4	100
2	0,1	0,9	50
3	0,1	1,9	25
4	0,1	3,9	12,5

- Mix without shaking.
- Patient Sample Preparation:
 - Prepare a 1+4 dilution of the patient plasma or control plasma in Owren's Buffer.
 - Mix without shaking.

- Testing:
 - Pipette, in duplicate, 0.1 mL of factor deficient plasma into a reaction tube.
 - Add 0.1 mL of standard, patient or control plasma dilution and incubate at *37°C for 2 minutes.
 - Add 0.1 mL of APTT Si L Minus reagent and incubate for 5 minutes at *37°C.
 - Add 0.1 mL of 0.025 M calcium chloride solution while simultaneously starting a stopwatch.
 - Determine the clot time for each of the standard, control or patient dilutions.
 - Plot % Activity (X-axis) versus Mean Clot Time (Y-axis) for the standards on 2 cycle log-log graph paper.
 - A straight line should be obtained.

Automated Method

Refer to the appropriate instrument operator manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application guides.

INTERPRETATION OF RESULTS

Interpolate patient values from the curve. Calculate exact patient or control values by correcting for differences in calibration plasma values as follows:

$$\text{Correct Activity (\%)} = \left(\frac{\text{Calibration plasma Reference Value}}{\text{Patient or Control Value}} \right) \times \text{Interpolated Patient or Control Value}$$

LIMITATIONS

The results obtained with Factor Deficient Plasmas depend on several factors strongly associated with instrumentation, types of reagents, deficient substrates and laboratory to laboratory variations^{3,4,5}. Each laboratory should establish an expected range for the particular instrument-reagent system.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. Helena Biosciences Europe supplies the following controls available for use with this product:

REF 5301	Speciality Assayed Control N (SAC-N)
REF 5302	Speciality Assayed Control A (SAC-A)

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range. Expected values for factor activity are 50–150%⁶.

PERFORMANCE CHARACTERISTICS

Each laboratory should establish its own performance data. Helena Biosciences Europe intrinsic factor assays are designed to give a linear standard curve from 10–150%. Within run and between run precisions are expected to be <5%, using a range of automated instruments.

BIBLIOGRAPHY

- Penner JA (1979) The University of Michigan Medical School Blood Coagulation Laboratory Manual, 14th Ed., University Publications, Ann Arbor, 72-78.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.

- Approved Guideline, 5th edn. CLSI: H21-A5.
- Kirkwood TBL *et al.* (1977) Identification of Sources of Variation in Factor VIII Assay, *British Journal of Haematology*, **37**:559-568.
- Goldenfarb MD (1971) Reproducibility in Coagulation Assays, *AJCP*, **55**:561-564.
- Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, *AJCP*, **59**:231-235.
- Triplett DA, Harms CS (1981) Procedures for the Coagulation Laboratory, Am. Society for Clin. Path, Chicago, 36

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BIBLIOGRAPHIE

- Penner JA (1979) The University of Michigan Medical School Blood Coagulation Laboratory Manual, 14th Ed., University Publications, Ann Arbor, 72-78.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
- Kirkwood TBL *et al.* (1977) Identification of Sources of Variation in Factor VIII Assay, *British Journal of Haematology*, **37**:559-568.
- Goldenfarb MD (1971) Reproducibility in Coagulation Assays, *AJCP*, **55**:561-564.
- Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, *AJCP*, **59**:231-235.
- Triplett DA, Harms CS (1981) Procedures for the Coagulation Laboratory, Am. Society for Clin. Path, Chicago, 36

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BIBLIOGRAFIA

- Penner JA (1979) The University of Michigan Medical School Blood Coagulation Laboratory Manual, 14th Ed., University Publications, Ann Arbor, 72-78.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
- Kirkwood TBL *et al.* (1977) Identification of Sources of Variation in Factor VIII Assay, *British Journal of Haematology*, **37**:559-568.
- Goldenfarb MD (1971) Reproducibility in Coagulation Assays, *AJCP*, **55**:561-564.
- Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, *AJCP*, **59**:231-235.
- Triplett DA, Harms CS (1981) Procedures for the Coagulation Laboratory, Am. Society for Clin. Path., Chicago, 36

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