Fibrinogen Calibrator

REF 5379

	Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom
	Tel: +44 (0)191 482 8440
	Fax: +44 (0)191 482 8442
	Email: info@helena-biosciences.com
	Web: www.helena-biosciences.com

HL-2-1390P 2016/01 (7)

Fibrinogen Calibrator Instructions for use

INTENDED PURPOSE

The Fibrinogen Calibrator kit is intended for use as a calibration material.

For use as a calibrator in the assay of fibrinogen in human plasma. May be used in conjunction with Helena Biosciences Europe Clauss Fibrinogen (Thrombin Only) (REF 5374, 5378) or Clauss Fibrinogen 100 (REF 5376, 5376H, 5376HR, 5376R).

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of:

Hepatitis B Antigen (HbsAg)
HIV 1 antibody
HIV 2 antibody
HCV antibody

However they should be handled with the same precautions as a human patient sample.

COMPOSITION

Component	Content	Description	Preparation
Fibrinogen Calibrator	10 x 1 mL	Each vial contains a lyophilised and stabilised preparation of pooled citrated human plasma.	Reconstitute each vial with 1 mL of distilled or deionised water, swirl gently and allow to stand for 10 minutes at room temperature. Mix well before use. Do not shake.

Each kit contains instructions for use.

Each kit contains lot specific reference values insert.

ITEMS REQUIRED BUT NOT PROVIDED

Fibrinogen Calibrator may be used when performing tests on any mechanical or photo-optical coagulation instrument in conjunction with all suitable, commercial reagents.

STORAGE, SHELF-LIFE AND STABILITY

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Fibrinogen Calibrator Avoid contamination. Once reconstituted, the plasma is stable for 4 hours at *2 –*8°C.

SAMPLE COLLECTION AND PREPARATION

Not applicable.

PROCEDURE

Each calibrator should be treated in the same manner as the unknown specimens and controls in accordance with the instructions outlined in each particular test protocol. Refer to the appropriate instrument operator manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application guides.

INTERPRETATION OF RESULTS

The assayed percent activities should be taken from the “Reference Value” column when Fibrinogen Calibrator is used to determine the laboratory’s standard curves. Ensure the lot number printed on this assay sheet is the same as that on the vial of Fibrinogen Calibrator to be used.

LIMITATIONS

The results obtained with Fibrinogen Calibrator depend on several factors strongly associated with instrumentation, types of reagents, deficient substrates and laboratory to laboratory variations^{1,2,3}. Each laboratory should establish an expected range for the particular instrument-reagent system.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges. Lot specific expected values are provided with each pack.

PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives using an opto-mechanical coagulation instrument. Each laboratory should establish its own performance data.

Reproducibility			
Inter-batch precision			
<i>Parameter</i>	<i>n</i>	<i>Mean</i>	<i>CV (%)</i>
Fibrinogen (g/L)	5	3.0	3.6

BIBLIOGRAPHY

- Kirkwood TBL *et al.* (1977) Identification of Sources of Variation in Factor VIII Assay, *British Journal of Haematology*, **37**:559-568.
- Goldenfarb MD (1971) Reproducibility in Coagulation Assays, *AJCP*, **55**:561-564.
- Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, *AJCP*, **59**:231-235.

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