ISI Calibrant Plasma Set

REF 5519

Helena Biosciences Europe Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 0SD, United Kingdom
Tel: +44 (0)191 482 8440
Fax: +44 (0)191 482 8442
Email: info@helena-biosciences.com
Web: www.helena-biosciences.com

HL-2-0446P 2015/10 (13)

ISI Calibrant Plasma Set Instructions for use

INTENDED PURPOSE

The ISI Calibrant Plasma Set kit is intended for use as a calibration material.

The ISI Calibrant Plasma Set is designed to simplify the laboratory calibration of thromboplastin International Sensitivity Index (ISI) values when using automated equipment for the determination of prothrombin time. The ISI Calibrant Plasma Set is intended for use with thromboplastin of human or rabbit origin and its performance has been validated with thromboplastin ISI's less than 2.00. The local ISI calibration of a thromboplastin enables laboratories to report International Normalised Ratio (INR) values much closer to the "true" INR value by correcting for coagulum effects on thromboplastin ISI. Other benefits of using the ISI Calibrant Plasma Set are a) they remove the requirement for running parallel PT times using the manual technique, b) a local supply of International Reference Preparation is not required and c) fresh plasmas from 60 long-term patients on therapeutic warfarin are not required. Thomson and colleagues also showed that both lyophilised plasmas from warfarin-treated patients and lyophilised plasmas artificially depleted warfarin-dependent clotting factors could be used to assess the effects of coagulometers on INR. Current guidelines based on a European Concerted Action on Anticoagulation (ECA) study showed that a minimum of 7 normal plasmas should be used for ISI calibration. The ISI Calibrant Plasma Set consists of 20 artificially depleted and 7 normal lyophilised plasmas that are calibrated against the WHO primary standards for Human recombinant and Rabbit derived thromboplastin. The ISI of a thromboplastin (TPT) which INR are derived is based upon a manual PT technique and can be considerably modified by coagulometers. The local system ISI of a thromboplastin can be calculated by comparing the PT times obtained for the ISI Calibrant Plasma Set plasmas on the thromboplastin / instrument system in laboratory use with the certified PT times obtained on the same lot of plasmas using the International Reference Preparation and the manual tit-tube technique.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 antibody HIV 2 antibody HCV antibody. However they should be handled with the same precautions as a human patient sample.

COMPOSITION

Component	Content	Description	Preparation
ISI Calibrant Plasma	27 x 0,5 mL	27 individually certified vials. Each vial contains either 0,5 mL of normal plasma or normal human plasma artificially depleted of clotting factors by barium sulphate adsorption with buffering agents (HEPES).	Reconstitute each plasma from the set with exactly 0,5 mL of distilled or deionised water. Allow to stand for 10 minutes and mix well before use. DO NOT SHAKE.

Each kit contains Instructions For Use.

Each kit contains Calibration Results Report Form.

ITEMS REQUIRED BUT NOT PROVIDED

- Any high quality electro-mechanical or photo-optical coagulation instrument suitable for the determination of prothrombin time (MPNT).
- Pipetting device to deliver 0,5 mL.
- 20 fresh plasma samples from normal individuals of both sexes with a wide age spread for the determination of the geometric mean value for the instrument used in conjunction with the ISI Calibrant Plasma Set.
- Thromboplastin Reagent, prepared according to the manufacturers' instructions.

STORAGE, SHELF-LIFE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Component	Stability	Signs of Deterioration
ISI Calibrant Plasma	Reconstituted plasmas are stable for 30 minutes at 15°-30°C or 4 hours at 2°-8°C. Storage longer than this time is not recommended.	The lyophilised plasmas should be dry, straw or pale-yellow coloured pellets or pieces of pellets. Failure to reconstitute to a homogeneous liquid, or the presence of particulate contamination may indicate product deterioration.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at 2°-8°C or 18°-24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at 37°C prior to testing. Do not keep at 37°C for more than 5 minutes.

PROCEDURE

IMPORTANT: DO NOT use the components of this kit with any other lot number kit. Each batch of ISI Calibrant Plasma Sets is individually certified.

- Reconstitute each plasma from the set as instructed in the "Composition" section.
- Perform four replicate prothrombin time tests on the vial using the appropriate thromboplastin and instrument combination. Refer to the appropriate instrument operator manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application notes.
- Record each individual PT value to the nearest 0.1 seconds together with the geometric mean value on the Calibration Sheet included with each kit. Ensure the full lot number of each plasma vial is recorded on the Calibration Sheet.
- Perform quadruplicate prothrombin time tests on a minimum of 20 fresh normal plasmas as instructed in the "Materials required but not provided" section and report the geometric mean PT value (MPNT) to the nearest 0.1 seconds on the Calibration Sheet. If a MPNT value has already been determined for the lot of thromboplastin and instrument in use, report this value instead.
- Complete the remaining details on the Calibration Sheet and fax or email the results to Helena Biosciences Europe using the number or address given.

INTERPRETATION OF RESULTS

The coefficient of variation (CV) between the 4 replicate PT times should be <5% for each plasma. The calibrated ISI value for the thromboplastin used on the instrument tested will be returned by fax or e-mail within 2 working days. The Local System ISI is the value of the slope of the line indicated by orthogonal regression analysis. Only calibrations where the CI slope is below the acceptable. If the CV of the calibration slope following orthogonal regression analysis is greater than 5%, the analysis should be repeated after the local MPNT has been checked. If the CV of the slope on the repeat analysis remains greater than 5%, the performance of the thromboplastin and instrument in use should be reviewed as these affect the precision of the calibration.

LIMITATIONS

It is important to calibrate human source thromboplastins with Human Reference Preparation values and rabbit source thromboplastins with Rabbit Reference Preparation values. Failure to ensure this may result in inaccurate ISI calibrations and unacceptable CV of slope values for the calibration.

QUALITY CONTROL

It is recommended to calibrate the local system ISI of thromboplastin using the ISI Calibrant Plasma Set at the introduction of a new thromboplastin lot number, following instrument service or if quality control values fall outside appropriate limits to ensure ongoing satisfactory system performance.

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges. This is particularly important for local ISI calibration.

PERFORMANCE CHARACTERISTICS

The reliability of the abnormal plasmas used in the manufacture of the ISI Calibrant Plasma Set has been evaluated in two major studies. 1) a national UK survey of users of the Manchester Reagent Thromboplastin¹ and 2) an international study². In the latter study, thirty-seven laboratories from the USA, Canada, France, Germany, Italy and Spain, utilising eleven types of commercial thromboplastins in 65 test systems, provided similar substantive evidence of the value of the calibrants in correcting for coagulum effects on INR across a series of instruments and reagents. In this study, the mean deviation from target INR for a series of 10 plasmas from long term warfarin patients was reduced from 14.41% (±15.28) to 1.04% (±7.51). The instrument thromboplastin combinations used in this study are summarised in the following table:

Instrument	Thromboplastin	Manufacturer ISI	Number of systems
KoaguLab	RecombiPlasTin	0.99	10
KoaguLab	Ortho Brain Thromboplastin	1.92	10
KoaguLab	Stago, France	1.79-2.04	2
KoaguLab	Thromboplastin C	2.80	1
MLA*	RecombiPlasTin	0.99	10
MLA*	Ortho Brain Thromboplastin	1.92	10
MLA*	Thromboplastin C	2.80	1
MLA*	Innovin	0.96-1.04	3
MLA*	IS, Baxter	#	2
MLA*	Thromboplastin C +	#	2
ACL**	RecombiPlasTin	0.99	13
ACL**	Ortho Brain Thromboplastin	1.92	13
ACL**	ACL	1.41-2.42	3
ACL**	Stago, France	1.79-2.04	1
ACL**	Simplastin XL	#	2
Coag-A-Mate	RecombiPlasTin	0.99	3
Coag-A-Mate	Ortho Brain Thromboplastin	1.92	3
Coag-A-Mate	Stago, France	1.79-2.04	2
Coag-A-Mate	Thrombomat	#	1
Behring BFA	RecombiPlasTin	0.99	1
Behring BFA	Ortho Brain Thromboplastin	1.92	1
Behring BFA	Thromborel S	#	1
* MLA = Medical Laboratory Automation			
** ACL = Automated Coagulation Laboratory			
# ISI ranged between 1.07-2.00			

BIBLIOGRAPHY

- Thomson JM, Taberner DA, Poller L (1990) Automation and prothrombin time: a United Kingdom field study of two widely used coagulometers. *J. Clin. Pathol.* 43:679-84
- Poller L, Barrowcliffe TW, van den Besselaar AMPH, Jespersen J, Tripodi A, Houghton D (1998) European Concerted Action on Anticoagulation. Minimum lyophilised plasma requirement for ISI calibration. *Amer. J. Clin. Pathol.* 109:196-204
- Clinical and Laboratory Standards Institute (2005) Procedures for Validation of INR and Local Calibration of PT/INR Systems. Approved Guidelines, CLSI H54-A, Vol. 25, No 23
- Poggio, M, van den Besselaar AMPH, van der Velde EA, Bertina RM (1998) The effect of some instruments for prothrombin time testing on the International Sensitivity Index (ISI) of two rabbit tissue thromboplastin reagents. *Thromb. Haemost.* 62: 868-74
- D'Angelo A, Seveso MP, D'Angelo SV, et al (1989) Comparison of two automated coagulometers and the manual tit tube method for the determination of prothrombin time. *Am. J. Clin. Pathol.* 92: 321-8
- Poller L, Thomson JM, Taberner DA (1989) Effect of automation on the prothrombin time test in NEQAS surveys. *J. Clin. Pathol.* 42: 97-100
- van Rijn JLM, Schmidt NA, Rutten W (1989) Correction of instrument and reagent based differences in determination of International Normalised Ratio (INR) for monitoring anticoagulant therapy. *Clin. Chem.* 35: 840-3

- Ray MJ, Smith IR (1990) The dependence of the International Sensitivity Index on the coagulometer used to perform the prothrombin time. *Thromb. Haemost.* 63: 424-9
- Clarke K, Taberner DA, Thomson JM, Morris JA, Poller L (1992) Assessment of value of calibrated lyophilised plasmas to determine International Sensitivity Index for coagulometers. *J. Clin. Pathol.* 45: 58-60
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Poller L, Thomson JM, Taberner DA, Clarke DK (1994) The correction of coagulometer effects on INR. A multicentre evaluation. *Brit. J. Haemat.* 86: 112-7
- Poller L, Houghton D, Carroll J (1994) The reliability of the MPNT of fresh plasmas and of the normal value from a lyophilised "normal" plasma in prothrombin time ratio determination. *Brit. J. Haemat.* 88: 866-73
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): an international multicentre study. *Amer. J. Clin. Pathol.* 103: 358-65
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios. *Amer. J. Clin. Pathol.* 103: 366-71

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BIBLIOGRAPHIE

- Thomson JM, Taberner DA, Poller L (1990) Automation and prothrombin time: a United Kingdom field study of two widely used coagulometers. *J. Clin. Pathol.* 43:679-84
- Poller L, Barrowcliffe TW, van den Besselaar AMPH, Jespersen J, Tripodi A, Houghton D (1998) European Concerted Action on Anticoagulation. Minimum lyophilised plasma requirement for ISI calibration. *Amer. J. Clin. Pathol.* 109:196-204
- Clinical and Laboratory Standards Institute (2005) Procedures for Validation of INR and Local Calibration of PT/INR Systems. Approved Guidelines, CLSI H54-A, Vol. 25, No 23
- Poggio, M, van den Besselaar AMPH, van der Velde EA, Bertina RM (1998) The effect of some instruments for prothrombin time testing on the International Sensitivity Index (ISI) of two rabbit tissue thromboplastin reagents. *Thromb. Haemost.* 62: 868-74
- D'Angelo A, Seveso MP, D'Angelo SV, et al (1989) Comparison of two automated coagulometers and the manual tit tube method for the determination of prothrombin time. *Am. J. Clin. Pathol.* 92: 321-8
- Poller L, Thomson JM, Taberner DA (1989) Effect of automation on the prothrombin time test in NEQAS surveys. *J. Clin. Pathol.* 42: 97-100
- van Rijn JLM, Schmidt NA, Rutten W (1989) Correction of instrument and reagent based differences in determination of International Normalised Ratio (INR) for monitoring anticoagulant therapy. *Clin. Chem.* 35: 840-3
- Ray MJ, Smith IR (1990) The dependence of the International Sensitivity Index on the coagulometer used to perform the prothrombin time. *Thromb. Haemost.* 63: 424-9
- Clarke K, Taberner DA, Thomson JM, Morris JA, Poller L (1992) Assessment of value of calibrated lyophilised plasmas to determine International Sensitivity Index for coagulometers. *J. Clin. Pathol.* 45: 58-60
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Poller L, Thomson JM, Taberner DA, Clarke DK (1994) The correction of coagulometer effects on INR. A multicentre evaluation. *Brit. J. Haemat.* 86: 112-7
- Poller L, Houghton D, Carroll J (1994) The reliability of the MPNT of fresh plasmas and of the normal value from a lyophilised "normal" plasma in prothrombin time ratio determination. *Brit. J. Haemat.* 88: 866-73
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): an international multicentre study. *Amer. J. Clin. Pathol.* 103: 358-65
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios. *Amer. J. Clin. Pathol.* 103: 366-71

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ISI variaiva tra 1,07 e 2,00

BIBLIOGRAFIA

- Thomson JM, Taberner DA, Poller L (1990) Automation and prothrombin time:a United Kingdom field study of two widely used coagulometers. *J. Clin. Pathol.*, 43:679-84
- Poller L, Barrowcliffe TW, van den Besselaar AMHP, Jespersen J, Tripodi A, Houghton D (1998) European Concerted Action on Anticoagulation. Minimum lyphised plasma requirement for ISI calibration. *Amer. J. Clin. Pathol.*, 109:196-204
- Clinical and Laboratory Standards Institute (2005) Procedures for Validation of INR and Local Calibration of PT/INR Systems: Approved Guidelines, CLSI H54-A, Vol. 25, No 23
- Poggio, M, van den Besselaar AMPH, van der Velde EA, Bertina RM (1988) The effect of some instruments for prothrombin time testing on the International Sensitivity Index (ISI) of two rabbit tissue thromboplastin reagents, *Thromb. Haemost.*, 62: 866-74
- D'Angelo A, Seveso MP, D'Angelo SV, et al (1989) Comparison of two automated coagulometers and the manual titr tube method for the determination of prothrombin time. *Am. J. Clin. Pathol.*, 92: 321-8
- Poller L, Thomson JM, Taberner DA (1989) Effect of automation on the prothrombin time test in NEQAS surveys, *J. Clin. Pathol.*, 42: 97-100
- van Rijn JLM, Schmidt NA, Rutten W (1989) Correction of instrument and reagent based differences in determination of International Normalised Ratio (INR) for monitoring anticoagulant therapy. *Clin. Chem.*, 355: 840-3
- Ray MJ, Smith IR (1990) The dependence of the International Sensitivity Index on the coagulometer used to perform the prothrombin time. *Thromb. Haemost.*, 63: 424-9
- Clarke K, Taberner DA, Thomson JM, Morris JA, Poller L (1992) Assessment of value of calibrated lyophilised plasmas to determine International Sensitivity Index for coagulometers. *J. Clin. Pathol.*, 45: 58-60
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Poller L, Thomson JM, Taberner DA, Clarke DK (1994) The correction of coagulometer effects on INR. A multicentre evaluation. *Brit. J. Haemat.*, 86: 112-7
- Poller L, Houghton D, Carroll J (1994) The reliability of the MNPT of fresh plasmas and of the normal value from a lyophilised "normal" plasma in prothrombin time ratio determination. *Brit. J. Haemat.*, 88: 866-73
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): an International multicentre study. *Amer. J. Clin. Pathol.*, 103: 358-65
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios. *Amer. J. Clin. Pathol.*, 103: 366-71

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* MLA = Medical Laboratory Automation

** ACL = Automated Coagulation Laboratory

EI ISI osciló entre 1,07-2,00

BIBLIOGRAFIA

- Thomson JM, Taberner DA, Poller L (1990) Automation and prothrombin time:a United Kingdom field study of two widely used coagulometers. *J. Clin. Pathol.*, 43:679-84
- Poller L, Barrowcliffe TW, van den Besselaar AMHP, Jespersen J, Tripodi A, Houghton D (1998) European Concerted Action on Anticoagulation. Minimum lyophilised plasma requirement for ISI calibration. *Amer. J. Clin. Pathol.*, 109:196-204
- Clinical and Laboratory Standards Institute (2005) Procedures for Validation of INR and Local Calibration of PT/INR Systems: Approved Guidelines, CLSI H54-A, Vol. 25, No 23
- Poggio, M, van den Besselaar AMPH, van der Velde EA, Bertina RM (1988) The effect of some instruments for prothrombin time testing on the International Sensitivity Index (ISI) of two rabbit tissue thromboplastin reagents, *Thromb. Haemost.*, 62: 866-74
- D'Angelo A, Seveso MP, D'Angelo SV, et al (1989) Comparison of two automated coagulometers and the manual titr tube method for the determination of prothrombin time. *Am. J. Clin. Pathol.*, 92: 321-8
- Poller L, Thomson JM, Taberner DA (1989) Effect of automation on the prothrombin time test in NEQAS surveys, *J. Clin. Pathol.*, 42: 97-100
- van Rijn JLM, Schmidt NA, Rutten W (1989) Correction of instrument and reagent based differences in determination of International Normalised Ratio (INR) for monitoring anticoagulant therapy. *Clin. Chem.*, 355: 840-3
- Ray MJ, Smith IR (1990) The dependence of the International Sensitivity Index on the coagulometer used to perform the prothrombin time. *Thromb. Haemost.*, 63: 424-9
- Clarke K, Taberner DA, Thomson JM, Morris JA, Poller L (1992) Assessment of value of calibrated lyophilised plasmas to determine International Sensitivity Index for coagulometers. *J. Clin. Pathol.*, 45: 58-60
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Poller L, Thomson JM, Taberner DA, Clarke DK (1994) The correction of coagulometer effects on INR. A multicentre evaluation. *Brit. J. Haemat.*, 86: 112-7
- Poller L, Houghton D, Carroll J (1994) The reliability of the MNPT of fresh plasmas and of the normal value from a lyophilised "normal" plasma in prothrombin time ratio determination. *Brit. J. Haemat.*, 88: 866-73
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): an International multicentre study. *Amer. J. Clin. Pathol.*, 103: 358-65
- Poller L, Triplett DA, Hirsh J, Carroll J, Clarke K (1995) A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios. *Amer. J. Clin. Pathol.*, 103: 366-71

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