LA Positive Control S

REF 5486
REF 5486D

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 0SD, United Kingdom
Tel: +44 (0)191 482 8440
Fax: +44 (0)191 482 8442
Email: info@helena-biosciences.com
Web: www.helena-biosciences.com

HL-2-0724P 2016/01 (9)

LA Positive Control S Instructions for use en

INTENDED PURPOSE

The LA Positive Control S kit is intended for use as a quality control material.

LA Positive Control S is prepared with plasma from human donors positive for Lupus anticoagulants. The plasma gives results typical of a Lupus Anticoagulant patient in DRVVT Screen, DRVVT Confirm and APTT-based tests.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg), HIV 1 antibody, HIV 2 antibody, HCV antibody. However they should be handled with the same precautions as a human patient sample.

COMPOSITION

Component	Content	Description	Preparation
LA Positive Control S	1 x 1 mL	Each kit contains lyophilised, citrated human plasma.	Reconstitute each vial with 1 mL of distilled/deionised water. Swirl gently and allow to stand for 15 minutes. Mix well before use. Do not shake.

Each kit contains Instructions for use.
Each kit contains lot specific reference values insert.

ITEMS REQUIRED BUT NOT PROVIDED

REF 5562 / 5558 / 5559 / 5562SLQ / 5558SLQ / 5559SLQ	APTT Si L Minus
REF 5560 / 5561	APTT Si L Plus
REF 5387ELQ / 5388ELQ	APTT SA
REF 5484	DRVVT Screen
REF 5485	DRVVT Confirm

STORAGE, SHELF-LIFE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label. Reconstituted vials of plasma should be kept on ice and are stable for 4 hours at $2-8^{\circ}\text{C}$.

SAMPLE COLLECTION AND PREPARATION

Not applicable

PROCEDURE

Each control should be treated in the same manner as the unknown specimen in accordance with the instructions outlined in each particular test protocol.

INTERPRETATION OF RESULTS

LA Positive Control S should show abnormal results when tested with DRVVT Screen reagent which are more normal when tested with the DRVVT Confirm reagent. Refer to the reference values insert for lot specific reference ranges.

LIMITATIONS

The results obtained with LA Positive Control S depends on several factors strongly associated with instrumentation, types of reagents, and laboratory to laboratory variations^{1,2,3}. Each laboratory should establish an expected range for the particular instrument-reagent system.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own range.

PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives as a guideline. Each laboratory should establish its own performance data.

Reproducibility								
PT-LI	APTT Si		DRVVT Screen		DRVVT Confirm			
<i>Clot formation (seconds)</i>	<i>CV (%)</i>	<i>Clot formation (seconds)</i>	<i>CV (%)</i>	<i>Clot formation (seconds)</i>	<i>CV (%)</i>	<i>Clot formation (seconds)</i>	<i>CV (%)</i>	<i>Ratio</i>
12,3	0,72	34,20	1,0	47,3	4,96	33,8	2,52	1,40

BIBLIOGRAPHY

- Kirkwood TBL *et al.* (1977) Identification of Sources of Variation in Factor VIII Assay, *British Journal of Haematology*, **37**:559-568.
- Goldenfarb MD (1971) Reproducibility in Coagulation Assays, *AJCP*, **55**:561-564.
- Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, *AJCP*, **59**:231-235.

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