# Manual D-Dimer

REF 5250

REF 5250H

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 0SD, United Kingdom

Tel: +44 (0)191 482 8440

Fax: +44 (0)191 482 8442

Email: info@helena-biosciences.com

Web: www.helena-biosciences.com

HL-2-0554P 2016/01 (14)

## Manual D-Dimer

Instructions for use

en

### INTENDED PURPOSE

The Manual D-Dimer kit is intended for carrying out immunoassay based haemostasis assays.

Manual D-Dimer is a latex agglutination test for semi-quantitative determination of D-dimer. D-dimer containing moieties are formed by fibrin degradation of Factor Xa cross-linked fibrin. Elevated levels of D-dimer are found in clinical conditions such as deep vein thrombosis (DVT), pulmonary embolism (PE) and other thromboembolic diseases. D-dimer levels are also elevated during pregnancy and high levels are associated with complications<sup>1</sup>. Manual D-Dimer utilizes a monoclonal antibody reacting with fibrin-D-dimer or fragment D of fibrin but not with intact fibrinogen<sup>2</sup>, allowing D-dimer determination in human plasma. Serum samples suited for FDP analysis can also be used. No international standard for D-dimer exists at present. For numerical values Helena has aligned the Manual D-Dimer with commercially available enzyme-linked immunosorbent assay (ELISA) tests.

### WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of:

Hepatitis B Antigen (HbsAg)

HIV 1 antibody

HIV 2 antibody

HCV antibody

However they should be handled with the same precautions as a human patient sample.

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Each kit contains instructions for use.

### ITEMS REQUIRED BUT NOT PROVIDED

- Pipettes for 200, 100 and 1000 µL volumes
- Pipette tips
- 3 mL test tubes
- 4 Timer

### STORAGE, SHELF-LIFE AND STABILITY

The reagents should be stored at  $2^{\circ}$ – $8^{\circ}$ C and used prior to the expiration date recorded on the box. Control plasmas, once reconstituted may be stored at  $2^{\circ}$ – $8^{\circ}$ C for one week or one month at  $-20^{\circ}$ C. The rubber stopper in the control vials may be discarded after visual inspection.

### SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes<sup>3</sup>. Plasma should be kept at  $2^{\circ}$ – $8^{\circ}$ C or  $18^{\circ}$ – $24^{\circ}$ C. Testing should be completed within 24 hours of sample collection, or plasma can be stored frozen at  $-20^{\circ}$ C or  $-70^{\circ}$ C for 24 months<sup>3</sup>. Thus quickly at  $37^{\circ}$ C prior to testing. Do not keep at  $37^{\circ}$ C for more than 5 minutes.

### PROCEDURE

Allow all vials to warm to room temperature before use, minimum 10 minutes.

#### A. Qualitative method

- Place 20 µL of sample, positive and negative control plasmas in circles on a test card.
- Place 20 µL of Manual D-Dimer latex suspension in a nearby area of each circle.
- Quickly mix the sample and Manual D-Dimer latex using clean mixing sticks for each sample. Start the timer.
- Rock the test card gently and read agglutination between 180 and 200 seconds. Positive (+) or negative (-) agglutination are compared to results obtained using the controls. The positive control is merely qualitative and should not be further diluted. Non-agglutinated latex means the sample is normal and no further testing is required.

#### B. Semi-quantitative method

Performed only on samples tested positive:

- Serially dilute 100 µL of sample 1:2, 1:4 and 1:8 with 100 µL Saline Solution using small test tubes.
- Mark the positions of sample dilutions on the test card and mix with latex suspension as per Qualitative method.
- The D-dimer concentration may be determined from the table in "Reference Values".

### INTERPRETATION OF RESULTS

1. Agglutination occurs within 180–200 seconds for samples containing more than 250 ng/mL D-dimer. The mean level of D-dimer in a healthy population is between about 8 and 135 ng/mL, and neat plasma from normal, healthy individuals should not agglutinate. The diagnostic predictive value of the Manual D-Dimer for thrombosis is high. The circulatory half-life of D-dimer is approximately 12 hours. Elevated D-dimer levels are therefore persistent for some time after the active process has ceased. In clinical studies on normal subjects, patients with phlebographically confirmed DVT, patients with DIC and patients with pre-eclampsia (Pre-EC) the following results were obtained; titers indicate the highest dilution at which the sample shows agglutination:

Sample	Number of patients	Negative D-Dimer	Plasma dilution	Undiluted	1:2	1:4	1:8
Normal	101	100		1*	-	-	-
DVT	48	3		10*	7*	14*	14*
DIC	29	0		3	3	4	19
Pre-EC	6	2		1	3	-	-

\* The agglutination was inhibited by addition of the D-dimer specific antibody (0.2 mg/mL) but not with a non-related antibody PAM-1.

2. The results may be reported either in D-dimer units or in fibrinogen equivalent units (FEU); 1 ng/mL of D-dimer is equivalent to approximately 2 ng/mL of FEU<sup>4</sup>.

3. Agglutination will be more pronounced and appears more rapidly at higher D-dimer concentrations.

Should you require quantitative D-Dimer analysis, Helena Biosciences Europe supplies the following kits for use on photo-optical instruments:

REF 5551 / 5552 Auto Blue D-Dimer 400

REF 5601 / 5501 Auto Red D-Dimer 700

### LIMITATIONS

1. A negative Manual D-Dimer test does not completely rule out thrombosis. A negative predictive value for patients with suspected DVT has been found to be 94% with the Manual D-Dimer kit<sup>5</sup>. Detection of elevated levels of D-dimer should be used with other clinical information in forming a diagnosis.

2. Agglutination with samples containing normal D-dimer levels may occasionally occur due to non-specificity.

3. A small number of samples, when mixed with the latex, may exhibit white flakes which should not be confused with agglutination.

### QUALITY CONTROL

The positive and negative controls provided should be used for quality control of each kit. It is recommended that both positive and negative controls are tested each time the kit is used. If either positive or negative control fails to elicit the appropriate response, patient results obtained on that occasion should not be used and the kit should be discarded or returned to Helena Biosciences Europe.

### REFERENCE VALUES

Agglutination occurs within 180 to 200 seconds for samples containing more than 250 ng/mL D-dimer. By testing serially diluted plasmas, semi-quantitative results can be obtained.

### D-Dimer Level (ng/mL)

ng/mL	Undiluted	1:2	1:4	1:8
<250	-	-	-	-
250-500	+	-	-	-
500-1000	+	+	-	-
1000-2000	+	+	+	-
>2000	+	+	+	+

Agglutination may be more pronounced and appears more rapidly at higher D-dimer concentrations. If you wish to express the results in fibrinogen equivalent units (FEU) multiply the D-dimer levels in the table above by 2, e.g. <250 ng/mL becomes <500 ng/mL. As discussed in several reports<sup>6,7</sup>, some commercial latex tests do not have the claimed sensitivity when compared to commercial ELISA.

### PERFORMANCE CHARACTERISTICS

#### Specificity

The monoclonal antibody used in this device is specific for D-dimer by virtue of the screening method used for hydridoma selection<sup>2</sup>. A hydridoma, secreting antibodies which reacted positively with purified D-dimer but not with whole fibrinogen or fragment D of fibrinogen, was selected. No cross-reactivity with fibrinogen or des-aa-fibrinogen was observed when either analysis was substituted for plasma in this assay. Plasma from 16 patients with rheumatoid arthritis was assayed and 14 of these were found to be non-agglutinating with the Manual D-Dimer test. The two agglutinations could be inhibited by addition of the D-dimer specific monoclonal antibody but not with addition of a monoclonal of the same subgroup, IgG1k (PAM-1). This suggests that the Manual D-Dimer is insensitive to rheumatoid factor disturbances.

#### Reproducibility

Three plasma samples were selected to test reproducibility of the assay. Each sample was tested 10 times on each of three different days. In each case, where a positive result was obtained, the sample was titrated. The results in D-dimer (ng/mL):

Sample	D-Dimer Level	Result
Normal	<250 ng/mL	Negative
Intermediate	3000 ng/mL	Titre 1:8
High	>16,000 ng/mL	Titre 1:64

Accuracy

The Manual D-Dimer kit was compared to another commercially available D-dimer latex. Both products gave a negative reaction when tested on 25 normal specimens. When 30 patient plasma samples were tested with ELISA and Manual D-Dimer, all samples with ELISA values above 225 ng/mL showed agglutination with Manual D-Dimer.

### BIBLIOGRAPHY

- Elms MJ *et al.* (1986) Rapid detection of cross-linked fibrin degradation products in plasma using monoclonal antibody-coated latex particles. *J. Clin. Path.* **49**:360-364
- Decker PV *et al.* (1987) Fibrinolytic response and fibrin fragment D dimer in patients with deep vein thrombosis. *Thromb. Haemostas.* **58**:1025-1029
- Ballegeer V *et al.* (1987) Binding properties of monoclonal antibodies against human fragment D-dimer of cross-linked fibrin to human plasma clots in an *in vivo* model in rabbits. *Thrombosis and Haemostasis*, **61**(2):307-313
- Holovat P *et al.* (1989) Binding properties of monoclonal antibodies against human fragment D-dimer of cross-linked fibrin to human plasma clots in an *in vivo* model in rabbits. *Thrombosis and Haemostasis*, **61**:1030-1032
- Ballegeer V *et al.* (1989) Binding properties of monoclonal antibodies against human fragment D-dimer of cross-linked fibrin to human plasma clots in an *in vivo* model in rabbits. *Thrombosis and Haemostasis*, **61**:1030-1032
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Clinical and Laboratory Standards Institute (2011) Quantitative D-dimer for the Exclusion of Venous Thromboembolic Disease: Approved Guideline. CLSI: H59-A
- Hansson PO *et al.* (1994) Can laboratory testing improve screening strategies for deep vein thrombosis at an emergency unit? *J. Intern. Med.* **235**:143-151
- Bounaoua H *et al.* (1989) Measurement of plasma D-dimer for diagnosis of deep venous thrombosis. *Am. J. Clin. Path.* **91**:82-85

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## BIBLIOGRAFIA

- Elms MJ *et al.* (1986) Rapid detection of cross-linked fibrin degradation products in plasma using monoclonal antibody-coated latex particles. *J. Clin. Path.* 39:360-364
- Declerck PV *et al.* (1987) Fibrinolytic response and fibrin fragment D dimer in patients with deep vein thrombosis. *Thromb. Haemostas.* 58:1025-1028
- Balegier V *et al.* (1987) Fibrinolytic response to venous occlusion and fibrin fragment D-dimer levels in normal and complicated pregnancy. *Thromb. Haemostas.* 58:1030-1032
- Holovet P *et al.* (1989) Binding properties of monoclonal antibodies against human fragment D-dimer of cross-linked fibrin to human plasma clots in an *in vivo* model in rabbits. *Thrombosis and Haemostasis.* 61(2):307-313
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H2-A5
- Clinical and Laboratory Standards Institute (2011) Quantitative D-dimer for the Exclusion of Venous Thromboembolic Disease: Approved Guideline, CLSI: H5-A
- Hansson PO *et al.* (1994) Can laboratory testing improve screening strategies for deep vein thrombosis at an emergency unit? *J. Intern. Med.* 236:143-151
- Bounaimeau H *et al.* (1989) Measurement of plasma D-dimer for diagnosis of deep venous thrombosis. *Am. J. Clin. Path.* 91:82-85

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