# Owren's Buffer

helena  
Biosciences Europe

REF 5375

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 0SD, United Kingdom

Tel: +44 (0)191 482 8440

Fax: +44 (0)191 482 8442

Email: info@helena-biosciences.com

Web: www.helena-biosciences.com

HL-2-0654P 2015/10 (6)

## Owren's Buffer Instructions for use

en

### INTENDED PURPOSE

The Owren's Buffer kit is intended for carrying out clot based haemostasis assays.

Owren's Buffer is a universal buffer used in quantitative coagulation studies on human plasma, typically Fibrinogen (Claus) and Factor Assays. Owren's Buffer is used to dilute reference plasma for the preparation of standard curves in Fibrinogen and Factor Assays and is also used as diluent for the patient plasmas for Fibrinogen or Factor quantitation. The buffer can be used on manual, semi-automated and fully automated methods for the performance of Fibrinogen and Factor Assays. The instructions supplied with the analyser and/or the test procedure should be followed for use.

### WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

### COMPOSITION

Component	Content	Description	Preparation
Owren's Buffer	10 x 25 mL	Each bottle contains an aqueous solution of 28.4 mM barbitol, 125.4 mM sodium chloride and 0.05% sodium azide as preservative pH 7.2-7.6.	The buffer is ready for use as packaged.

Each kit contains Instructions For Use.

### ITEMS REQUIRED BUT NOT PROVIDED

Helena Biosciences Europe supply the following kits available for use with this product:

REF 5376/35	Clauss Fibrinogen 35
REF 5556	Clauss Fibrinogen 50
REF 5376	Clauss Fibrinogen 100

### STORAGE, SHELF-LIFE AND STABILITY

Unopened bottles are stable until the given expiry date when stored under conditions indicated on the bottle label. Avoid contamination of the buffer and discard if any particulate matter indicative of contamination is observed.

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