Protein C (Chromogenic)

REF 5543

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom
 Tel: +44 (0)191 482 8440
 Fax: +44 (0)191 482 8442
 Email: info@helena-biosciences.com
 Web: www.helena-biosciences.com

HL-2-1677P 2016/01 (7)

Protein C (Chromogenic) Instructions for use

INTENDED PURPOSE

The Protein C (Chromogenic) kit is intended for carrying out chromogenic based haemostasis assays.

Protein C is a vitamin K dependent protein which plays an important role in the regulation of anticoagulant mechanisms. It can be synthesized by inactivating factors Va and VIII or when activated, can stimulate fibrinolysis. Protein C circulates as a zymogen, and is converted to an active serine protease by the action of thrombin in the presence of thrombomodulin. Both hereditary and acquired Protein C deficiencies have been shown to be a risk factor for development of venous thrombosis. Protein C in plasma is activated by a specific fraction from the *Agkistrodon contortrix* snake venom. The amount of activated protein C (APC) is determined by monitoring the rate of hydrolysis of a protein C specific chromogenic substrate. The release of pNA is measured at 405 nm and is proportional to the protein C level. The Helena Biosciences Europe Protein C (Chromogenic) assay is intended for the quantitative determination of protein C in human plasma using a chromogenic assay method.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product label for appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

COMPOSITION

Component	Content	Description	Preparation
Protein C Substrate	6 x 2 mL	Each vial contains 2.75 µmol lyophilised pyro-Glu-Pro-Arg-pNA.HCl.	Reconstitute each vial with 2 mL of Protein C Diluent. If cloudy, warm at 37°C for a few minutes.
Protein C Activator	6 x 2 mL	Each vial contains 0.8 units of activator from snake venom (Protac [®]).	Reconstitute each vial with 2 mL deionised water.
Protein C Diluent	3 x 5 mL	Each vial contains tris buffer with sodium Azide as a preservative.	Ready for use.

ITEMS REQUIRED BUT NOT PROVIDED

- Coagulation instrument or spectrophotometer operable at 405 nm.
- REF 5185 Calibration Plasma.
- Glacial acetic acid.
- 37°C water bath or dry bath.
- Laboratory timer

STORAGE, SHELF-LIFE AND STABILITY

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Protein C Substrate	Reconstituted reagent is stable for 1 week at 2–8°C or one month at -20°C.
Protein C Activator	Reconstituted reagent is stable 1 week at 2–8°C or one month at -20°C.
Protein C Diluent	Store at 2–8°C.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at 2–8°C or 15–25°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at 37°C prior to testing. Do not keep at 37°C for more than 5 minutes. Erroneous results may be caused by contamination with tissue fluids or stasis. Avoid agitation, air bubbles or foaming. For the effects of commonly administered drugs, refer to Young, *et al.*

PROCEDURE

Prepare plasma standards (REF 5185) and patients' plasma samples as follows:

Important: Use only 0.85% NaCl solutions for dilutions.

Standard %	Plasma	Buffer
100%	100 µL Calibration Plasma	+ 300 µL saline
50%	50 µL Calibration Plasma	+ 350 µL saline
0%	400 µL saline only	+ 300 µL saline
0%	100 µL plasma	+ 300 µL saline

A. End Point Method Assay

- To a glass or plastic test tube:
 - Add 100 µL standard or patients' plasma dilution.
 - Incubate at 37°C for 2 minutes.
 - Add 200 µL Protein C Activator and mix.
 - Incubate at 37°C for 5 minutes.
 - Add 200 µL Protein C Substrate and mix.
 - Incubate at 37°C for 10 minutes
 - Add 200 µL acetic acid and mix.
 - * Add 200 µL water (optional).

* Some spectrophotometers require a minimum of 1 mL volume in the cuvette.

- Read absorbance at 405 nm in a 1 cm semi-micro cuvette against a blank prepared with deionised water. As many as ten determinations can be performed simultaneously with the same spot watch by staggering pipetting steps at five second intervals.


B. Kinetic Method

Automated analyser required. Refer to the appropriate Instrument Operator Manual for detailed instructions or contact Helena Biosciences Europe for instrument-specific application guides. Example table:

Reagent	Volume	Time
Sample	75 µL	3:00
Protein C Activator	75 µL	5:00
Protein C Substrate	75 µL	5:00

Calibration Curve

Plot the absorbance obtained for each of the protein C standards against protein C % on linear graph paper. The protein C concentration in patients' plasma specimens can be determined by interpolation from the calibration curve. If a commercial protein C standard is used, the protein C concentration in the patient's specimen should be adjusted for the protein C concentration in the standard. The calibration curve shown below is an example only. A calibration curve must be obtained each time the assay is performed.

INTERPRETATION OF RESULTS

Protein C deficiencies, either congenital or acquired, may lead to serious thrombotic events such as thrombophlebitis, deep vein thrombosis, or pulmonary embolism. Approximately 2-8% of all patients with venous thrombosis under the age of 40-45 years have Protein C deficiencies. Patients with hereditary deficiencies generally present with venous thrombosis in young adulthood. The first episode is usually spontaneous and associated with trauma or stress to the haemostasis mechanism. The prevalence of Protein C deficiencies is one case per 300 people or approximately 0.33%.

Acquired Deficiencies

Decreased levels of protein C are observed in the following cases:

- Hepatic disorders: hepatitis, cirrhosis.
- DIC.
- Oral anticoagulant therapies, in which cases the interpretation of test results is difficult if the patients have had a history of thromboses and are receiving oral anticoagulant treatments.

Congenital Deficiencies

Two types of deficiencies can be observed:

- Quantitative type I, in which the functional and antigenic levels are simultaneously decreased.
- Qualitative or type II, characterised by a decreased functional level and by a normal or subnormal immunological level.

LIMITATIONS

Icteric, haemolysed or hyperlipemic specimens interfere with absorbance readings, thus requiring the use of plasma blanks for accurate results. Plasma blanks are also needed on patients when contact factor activation is suspected, such as DIC patients, or from individuals on oral contraceptives where clot activation may occur. Blanks are obtained by substituting saline for the Protein C Activator in the test reaction. Subtract the sample blank activity from the sample test activity. The naturally-occurring activator of protein C is thrombin in the presence of thrombomodulin. The possible existence of clinical protein C abnormalities not detectable by reactivity with the snake venom activator or with the chromogenic substrate used in this kit should not be precluded. To assure accurate, reproducible results, use accurate pipetting devices and observe recommended procedures with particular emphasis on incubation times and incubation temperature.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, Helena Biosciences Europe supply the following controls available for use with this product:

REF 5301	Speciality Assayed Control N
REF 5302	Speciality Assayed Control A

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range. Protein C activity values are usually expressed in relative percentages compared to a pooled normal plasma standard. Bertina *et al.* reported a range of 65-145% in healthy individuals with diminished levels found following anticoagulant therapy. There is apparently no difference in protein C between healthy males and females. The following reference values have been determined by Helena Biosciences Europe or their representatives using a photo-optical coagulation instrument, using 17 plasmas of presumed healthy men and women.

Protein C activity %	SD	Range
105%	26.6	61.2-162.9%

PERFORMANCE CHARACTERISTICS

Helena Biosciences Europe or their representatives have determined the following performance characteristics as a guideline. Each laboratory should establish its own performance data.

Specificity

In studies where recovery of protein C plasma was determined after addition of various amounts of protein C to protein C deficient plasma, the expected recovery was obtained.

Reproducibility

n	Intra-assay precision		Inter-assay precision	
	Protein C	CV (%)	n	Protein C CV (%)
10	68.5	2.5	18	45.3 2.6

BIBLIOGRAPHY

- Walker FJ (1981) Regulation of activated protein C by protein S. The role of phospholipid in factor Va inactivation, *J Biol Chem*, **256**: 11128-31
- Taylor FB, Lockhart MS (1985) A new function for activated protein C: activated protein C prevents inhibition of plasminogen activators by release from mononuclear leucocytes-platelet suspensions stimulated by phorbol diester, *Thromb Res*, **37**(1): 155-64
- Fulcher CA, Gardiner JE, Griffin JH, Zimmerman TS (1984) Proteolytic inactivation of human factor VIII procoagulant protein by activated human protein C and its analogy with factor V, **63**(2): 486-9
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Young DS *et al.* Effects of Drugs on Clinical Laboratory Tests, 3rd ed. AC&C Press, Washington, D.C., 1990
- Bertina RM, Broekmans AW, Krommenhoek-Van Es C, Van Wijngaarden A (1984) The Use of a Functional and Immunologic Assay for Plasma Protein C in the Study of the Heterogeneity of Congenital Protein C Deficiency, *Thromb Haemostas*, **51**: 1-5
- Pabinger-Fasching I, Bertina RM, Lechner K, Nieser H, Korninger C (1983) Protein C Deficiency in Two Austrian Families, *Throm Haemostas*, **50**: 810-3

Protein C (Chromogenic) Fiche technique

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BIBLIOGRAFIA

- Walker FJ (1981) Regulation of activated protein C by protein S. The role of phospholipid in factor Va inactivation, *J Biol Chem*, **256**: 11128-31
- Taylor FB, Lockhart MS (1985) A new function for activated protein C: activated Protein C prevents inhibition of plasminogen activators by releaseate from mononuclear leukocytes–platelet suspensions stimulated by phorbol diester, *Thromb Res*, **37**(1): 155-64
- Fulcher CA, Gardiner JE, Griffin JH, Zimmerman TS (1984) Proteolytic inactivation of human factor VIII procoagulant protein by activated human protein C and its analogy with factor V, **63**(2): 486-9
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Young DS *et al.* Effects of Drugs on Clinical Laboratory Tests, 3rd ed., AACCP Press, Washington, D.C., 1990
- Bertina RM, Broekmans AW, Krommenhoek-Van Es C, Van Wijngaarden A (1984) The Use of a Functional and Immunologic Assay for Plasma Protein C in the Study of the Heterogeneity of Congenital Protein C Deficiency, *Thromb Haemostas*, **51**: 1-5
- Pabinger-Fasching I, Bertina RM, Lechner K, Niesser H, Korninger C (1983) Protein C Deficiency in Two Austrian Families, *Throm Haemostas*, **50**: 810-3

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