Ristocetin Cofactor Kit

REF 5370

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom

Tel: +44 (0)191 482 8440

Fax: +44 (0)191 482 8442

Email: info@helena-biosciences.com

Web: www.helena-biosciences.com

HL-2-0593P 2015/10 (11)

Instructions for use

INTENDED PURPOSE

The Ristocetin Cofactor Kit is intended for the determination of the % activity of von Willebrand Factor.

It is used to measure the ability of a patient's plasma to agglutinate formalin-fixed platelets in the presence of Ristocetin. The rate of ristocetin-induced platelet agglutination is related to the concentration of von Willebrand Factor. Helena Biosciences Europe's Ristocetin Cofactor Kit is intended for use in the quantitative determination of von Willebrand Factor activity in plasma .

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of: Hepatitis B Antigen (HbsAg) HIV 1 antibody HIV 2 antibody HCV antibody However they should be handled with the same precautions as a human patient sample.

COMPOSITION

Allow all vials to warm to room temperature before reconstitution.

Component	Content	Description	Preparation
Lyophilised Platelets	4 x 5 mL	Each vial contains lyophilised formalin fixed platelets and stabilisers.	Reconstitue each vial with 5 mL of Tris-Buffered Saline. Swirl gently and allow to mix for 20 minutes on a roller mixer. Mix well before each use. Do not shake.
Ristocetin: 10 mg/mL	2 x 1,5 mL	Contains 10 mg/mL Ristocetin lyophilised with stabilisers.	Reconstitute each vial with 1.5 mL of purified water. Swirl gently and allow to stand for 10 minutes. Mix gently before use. Do not shake.
Calibration Plasma	2 x 1 mL	Contains 1 mL of lyophilised normal human plasma assayed for Ristocetin Cofactor activity.	Reconstitute each vial with exactly 1 mL of purified water. Swirl gently and allow to stand for 10 minutes. Mix gently before use. Do not shake.
Ristocetin Cofactor Abnormal Control	2 x 0,5 mL	Contains 0,5 mL of lyophilised plasma from an individual with von Willebrand disease. The activity is printed on an insert sheet provided in the kit.	Reconstitute each vial with 0,5 mL of purified water. Swirl gently and allow to stand for 10 minutes. Mix gently before use. Do not shake.
Tris-Buffered Saline	1 x 35 mL	Contains 35 mL of buffer containing 0,05 M Tris and 0,15 M sodium chloride with sodium azide as preservative.	The buffer is ready for use as packaged.

Each kit contains instructions for use.

Each kit contains lot specific reference values insert, and a report form.

ITEMS REQUIRED BUT NOT PROVIDED

Platelet Aggregometer.

STORAGE, SHELF-LIFE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Lyophilised Platelets	Once reconstituted, the reagent is stable for 1 day (24 hours) at *2 –*8°C.
Ristocetin: 10 mg/mL	Once reconstituted, the reagent is stable for 8 hours at *2 –*8°C or 30 days at -20°C.
Calibration Plasma	Once reconstituted, the reagent is stable for 2 hours at *2 –*8°C or 30 days at -20°C.
Ristocetin Cofactor Abnormal Control	Once reconstituted, the reagent is stable for 8 hours at *2 –*8°C or 30 days at -20°C.
Tris-Buffered Saline	Opened bottle must be stored at *2 –*8°C.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at *18 –*24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at *37°C prior to testing. Do not keep at 37°C for more than 5 minutes⁵.

PROCEDURE

Prepare all reagents as instructed with each kit. Set up platelet aggregometer as indicated in the appropriate instruction manual. For Helena AggRAM please refer to the AggRAM step by step guide HL-2-1785P 2008-01(1). **A basic aggregometer protocol is outlined below:**

- Standard Curve Preparation:
 - Prepare the following dilutions in Tris-buffered saline:

Tube	Calibration Plasma (μL)	Tris-Buffered Saline (μL)	Activity (%)
1	100	100	100
2	100	300	50
3	100	700	25
4	100	1500	12,5
 - Mix dilutions thoroughly.
- Patient Sample Preparation:
 - Prepare 1+1 and 1+3 dilutions of the patient plasma or control plasma in Tris-buffered saline. Mix without shaking.
- Aggregation Blank Preparation:
 - Pipette 250 μL of platelet suspension and 250 μL of Tris-buffered saline into an aggregometer cuvette with a stir-bar. This is used to set the 100% activity level.
- Testing:
 - Set the 100% activity in each channel using the Aggregation Blank sample.
 - Pipette 400 μL of platelet suspension into an aggregometer cuvette for each sample to be tested. Add a stir-bar to each cuvette.
 - Place the cuvette into the aggregometer and add 50 μL of Ristocetin to the cuvette. Incubate at 37°C for exactly 3 minutes.
 - Add 50 μL of standard, patient or control plasma dilution.
 - Monitor platelet aggregation for 5 minutes.

INTERPRETATION OF RESULTS

Manual method

Plot Calibration Plasma dilution activity (X-axis) versus Aggregation Slope (Y-axis) on log-log graph paper. Draw a straight line of best fit. Interpolate patient slope values from the graph to give % activity. Correct these results for dilution and assayed activity of the Calibration Plasma.

Automated method

Appropriate aggregometer software should automatically plot standard curve and report results as final activity %.

LIMITATIONS

Avoid icteric, lipaemic, and haemolysed samples. The Ristocetin Cofactor Kit fails to reflect accurately von Willebrands' disease in several situations such as pregnancy⁶, infusion of commercial factor VIII concentrates^{7,8} or administration of DDaVP⁹. In such instances, vWF R:Co may be corrected, yet the bleeding time remains prolonged. In addition, vWF R:Co levels may be normal in Type II B von Willebrands' disease even though the bleeding time is prolonged⁷.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. Helena Biosciences Europe supplies the following controls available for use with this product:

REF 5301	Speciality Assayed Control N
REF 5302	Speciality Assayed Control A
REF 5373	Ristocetin Cofactor Abnormal Control

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges. Example values for Ristocetin Cofactor Activity in adults using the Helena AggRAM are 69–116% (0.69–1.16 units/mL).

PERFORMANCE CHARACTERISTICS

Performance conforms to CLSI H51-A, Assays of von Willebrand Factor Antigen and Ristocetin Cofactor; Approved Guideline, with Coefficients of Variation (CV) of 15% or less⁵.

BIBLIOGRAPHY

- Zimmerman TS, Ruggeri ZM, 'von Willebrands' Disease', Progress in Thrombosis and Haemostasis, Spaett (Ed), 1982, 6: 203-236.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
- Ratnoff OD, Bennett B (1973) Clues to Pathogenesis of Bleeding in von Willebrands' Disease, *N. Eng. J. Med.* **289**:1182-1183.
- Blatt PM *et al.* (1976) Antithemophilic factor Concentrate Therapy in von Willebrands' Disease, *J. Am. Med. Assn.* **236**:2770-2772.
- Green D., Potter EV (1976) Failure of AHF Concentrate to Control Bleeding in von Willebrands' Disease, *Am. J. Med.* **60**:357-360.
- Mannucci PM *et al.* (1976) Studies on the Prolonged Bleeding Time in von Willebrands' Disease, *J. Lab. Clin. Med.* **88**:662-671.
- Ruggeri ZM *et al.* (1980) Heightened Interaction Between Platelets and Factor VIII/von Willebrand Factor in a New Subtype of von Willebrands' Disease, *N. Eng. J. Med.* **302**:1047-1051
- Assays of von Willebrand Factor Antigen and Ristocetin Cofactor Activity; Approved Guideline, H51-A; Vol. 22, No. 20.

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BIBLIOGRAPHIE

- Zimmerman TS, Ruggeri ZM, 'von Willebrands' Disease', Progress in Thrombosis and Haemostasis, Spaett (Ed), 1982, 6: 203-236.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
- Ratnoff OD, Bennett B (1973) Clues to Pathogenesis of Bleeding in von Willebrands' Disease, *N. Eng. J. Med.* **289**:1182-1183.
- Blatt PM *et al.* (1976) Antithemophilic factor Concentrate Therapy in von Willebrands' Disease, *J. Am. Med. Assn.* **236**:2770-2772.
- Green D., Potter EV (1976) Failure of AHF Concentrate to Control Bleeding in von Willebrands' Disease, *Am. J. Med.* **60**:357-360.
- Mannucci PM *et al.* (1976) Studies on the Prolonged Bleeding Time in von Willebrands' Disease, *J. Lab. Clin. Med.* **88**:662-671.
- Ruggeri ZM *et al.* (1980) Heightened Interaction Between Platelets and Factor VIII/von Willebrand Factor in a New Subtype of von Willebrands' Disease, *N. Eng. J. Med.* **302**:1047-1051.
- Assays of von Willebrand Factor Antigen and Ristocetin Cofactor Activity; Approved Guideline, H51-A; Vol. 22, No. 20.

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BIBLIOGRAFIA

- Zimmerman TS, Ruggeri ZM, von Willebrands' Disease¹, Progress in Thrombosis and Haemostasis, Spaett (Ed), 1982, 6: 203-236.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5.
- Ratnoff OD, Bennett B (1973) Clues to Pathogenesis of Bleeding in von Willebrands' Disease, *N. Eng. J. Med.* **289**:1182-1183.
- Blatt PM *et al.* (1976) Antithemophilic factor Concentrate Therapy in von Willebrands' Disease, *J. Am. Med. Assn.* **236**:2770-2772.
- Green D., Potter EV (1976) Failure of AHF Concentrate to Control Bleeding in von Willebrands' Disease, *Am. J. Med.* **60**:357-360.
- Mannucci PM *et al.* (1976) Studies on the Prolonged Bleeding Time in von Willebrands' Disease, *J. Lab. Clin. Med.* **88**:662-671.
- Ruggeri ZM *et al.* (1980) Heightened Interaction Between Platelets and Factor VIII/von Willebrand Factor in a New Subtype of von Willebrands' Disease, *N. Eng. J. Med.* **302**:1047-1051.
- Assays of von Willebrand Factor Antigen and Ristocetin Cofactor Activity; Approved Guideline, H51-A; Vol. 22, No. 20.

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