# Coagulation Control Plasmas

**Helena**  
Biosciences Europe

REF 5186  
REF 5187  
REF 5183  
REF 5482

Routine Control N  
Routine Control A  
Routine Control SA  
Routine Coagulation Control Set

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 0SD, United Kingdom  
Tel: +44 (0)191 482 8440  
Fax: +44 (0)191 482 8442  
Email: info@helena-biosciences.com  
Web: www.helena-biosciences.com

HL-2-0482P 2016/01 (16)

## Coagulation Control Plasmas

en

### INTENDED PURPOSE

The Coagulation Control Plasma kit is intended for use as a quality control material. Routine Control N, Routine Control A and Routine Control SA are for use as normal, moderately prolonged and markedly prolonged controls for PT and aPTT assays. They are also assayed for Fibrinogen, TCT and ATIII, and are prepared from normal human plasma.

### WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations. Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of: Hepatitis B Antigen (HbsAg) Hepatitis C Antigen (HcAg) HIV 2 antibody However they should be handled with the same precautions as a human patient sample.

### COMPOSITION

REF	Component	Content	Description
5186	Routine Control - N	10 x 1 mL	Prepared from pooled normal plasma.
5187	Routine Control - A	10 x 1 mL	Prepared from asclerotic human plasma.
5183	Routine Control - SA	10 x 1 mL	Prepared from asclerotic human plasma.
5482	Routine Coagulation Control Set		
	Routine Control - N	4 x 1 mL	
	Routine Control - A	3 x 1 mL	
	Routine Control - SA	3 x 1 mL	

Each kit contains instructions for use.

Each vial contains 1 mL of buffered, lyophilised human plasma. Reconstitute each vial of the appropriate control with 1 mL of distilled or deionised water. Swirl gently. Allow to stand for 10 minutes for complete dissolution and mix well before use.

### ITEMS REQUIRED BUT NOT PROVIDED

Coagulation Control Plasmas may be used when performing tests on any mechanical or photo-optical coagulation instrument in conjunction with suitable commercial reagents.

### STORAGE SHELF-LIFE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label. The reconstituted controls are stable for 8 hours when kept at 2 – 8°C or 4 weeks at -20°C when flash frozen. Keep covered.

### SAMPLE COLLECTION AND PREPARATION

Not applicable.

### PROCEDURE

Each control should be treated in the same manner as the unknown specimen in accordance with the instructions outlined in each particular test protocol.

### INTERPRETATION OF RESULTS

Routine Control N should give values within the laboratory normal range for PT, aPTT and fibrinogen assays. Routine Control A and Routine Control SA have been standardised to give prolonged and markedly prolonged PT and aPTT times respectively. Lot and routine specific expected values are provided with each pack of controls.

### LIMITATIONS

The results obtained with Coagulation Control Plasmas depend on several factors strongly associated with instrumentation. Types of reagent, reagent substrate and laboratory to laboratory variations. Each laboratory should establish an expected range for the particular instrument being used.

### QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

### REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

### PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives using an opto-mechanical coagulation instrument. Each laboratory should establish its own performance data.

Reproducibility	Intra-assay precision	PT CV (%)	aPTT CV (%)
Sample	n		
Routine Control N	5	2/83	1/01
Routine Control A	5	2/76	1/71
Routine Control SA	5	1/72	1/03

### BIBLIOGRAPHY

1. Kirkwood TBL, et al. (1977) Identification of Sources of Variation in Factor VIII Assay. *British Journal of Haematology*, 37:555-568
2. Goldfarb MD (1971) Reproducibility in Coagulation Assays. *AJCP* 55:561-564.
3. Pallitt HA and Longbery JR (1979) A Precision Study of Coagulation Factor Assay Techniques. *AJCP* 59:231-235.

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