Speciality Assayed Controls

REF 5301 Speciality Assayed Control N
REF 5302 Speciality Assayed Control A

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom

Tel: +44 (0)191 482 8440

Fax: +44 (0)191 482 8442

Email: info@helena-biosciences.com

Web: www.helena-biosciences.com

HL-2-0427P 2016/01 (15)

 	
	Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom
	Tel: +44 (0)191 482 8440
	Fax: +44 (0)191 482 8442
	Email: info@helena-biosciences.com
	Web: www.helena-biosciences.com

HL-2-0427P 2016/01 (15)

	
--	--

Speciality Assayed Controls

INTENDED PURPOSE

The Speciality Assayed Controls kit is intended for use as a quality control material.

Speciality Assayed Control N (SAC - N) may be used as a normal control, and Speciality Assayed Control A (SAC - A) may be used as an abnormal control when assaying for Factors II, V, VII, VIII, IX, X¹, XI, XII², fibrinogen³, von Willebrand Factor, antigenic and functional Protein C and Protein S (total and free), as well as the chromogenic assays including Antithrombin III, Protein C, Factor VIII and Plasminogen. Speciality Assayed Controls are used in factor assays and other testing in the same manner as a fresh plasma pool. Factor II, VII, VIII, IX and X values and the chromogenic Factor VIII, Antithrombin III and Protein C values are traceable to World Health Organisation standards to ensure the utmost credibility in values stated⁴. The control plasma should be used to gauge internal factors in each laboratory's system.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of:
Hepatitis B Antigen (HbsAg)
HIV 1 antibody
HIV 2 antibody
HCV antibody
However they should be handled with the same precautions as a human patient sample.

REF	Component	Content	Description	Preparation
5301	SAC - N	10 x 1 mL	Prepared from a frozen pool of citrated plasma and is buffered and lyophilised to ensure stability of all plasma constituents ⁵ .	Reconstitute with 1 mL of distilled or deionised water. Swirl gently. Allow product 20 minutes for complete dissolution.
5302	SAC - A	10 x 1 mL	Prepared from a frozen pool of citrated plasma and is buffered and lyophilised to ensure stability of all plasma constituents ⁵ .	Reconstitute with 1 mL of distilled or deionised water. Swirl gently. Allow product 20 minutes for complete dissolution.

Each kit contains Instructions For Use.

Each kit contains lot specific reference values insert.

ITEMS REQUIRED BUT NOT PROVIDED

Speciality Assayed Controls may be used when performing tests on any mechanical or photo-optical coagulation instrument in conjunction with all suitable, commercial reagents.

STORAGE, SHELF-LIFE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Component	Stability	Signs of Deterioration
SAC - N (REF 5301) and SAC - A (REF 5302)	Values for Factor VIII, von Willebrand factor and ristocetin co-factor are stable for 2 hours at *2 –*8°C. All other factors are stable for 4 hours at *2 –*8°C or 4 weeks at -20°C when flash frozen.	Unreconstituted Speciality Assayed Controls must appear as a light yellow, dry plug. If any unusual conditions are noted, the customer should notify Helena Biosciences Europe before using.

SAMPLE COLLECTION AND PREPARATION

Not applicable.

PROCEDURE

Each control should be treated in the same manner as the unknown specimen in accordance with the instructions outlined in each particular test protocol.

INTERPRETATION OF RESULTS

The assayed percent activities of the various coagulation factors should be used as a guideline. If the product does not perform as expected, the patient results should be considered invalid. Ensure the lot number printed on this assay sheet is the same as that on the vial of Speciality Assayed Control to be used.

LIMITATIONS

The results obtained with Speciality Assayed Controls depend on several factors strongly associated with instrumentation, types of reagents, deficient substrates and laboratory to laboratory variations^{6,7,8}. Each laboratory should establish an expected range for the particular instrument-reagent system used.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

REFERENCE VALUES

Refer to the reference values insert for lot specific reference ranges. Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives using an opto-mechanical coagulation instrument. Each laboratory should establish its own performance data.

Inter-assay precision		SAC - N		SAC - A	
<i>Parameter</i>	<i>n</i>	<i>Mean</i>	<i>CV (%)</i>	<i>Mean</i>	<i>CV (%)</i>
Fibrinogen (g/L)	3	2,9	6,9	1,24	3,35
Factor IX (%)	3	132,7	4,3	50,3	5,4
Protein S (%)	3	110,4	5,0	42,4	3,1

BIBLIOGRAPHY

- Babson AL and Flanagan ML (1975) Quantitative One Stage Assays for Factors V and X, *AJCP*, **64**:817-819
- Hardisty RM, *et al.* (1962) A One Stage Factor VIII Assay and Its Use on Venous and Capillary Plasma, *Thrombosis et Diathesis Haemorrhagica*, **7**:215-229
- Morse EE, *et al.* (1971) Automated Fibrinogen Determination, *AJCP*, **55**:671-676
- Elodi S, Katalin V, Hollan S (1978) Some Sources of Error in the One-Stage Assay of Factor VIII, *Haemostasis*, **7**:1-9
- Thelin M (1968) Preparation and Standardization of a Stable AHF Plasm, *Thrombosis et Diathesis Haemorrhagica*, **19**:423
- Kirkwood TBL et al. (1977) Identification of Sources of Variation in Factor VIII Assay, *British Journal of Haematology*, **37**:559-568
- Goldenfarb MD (1971) Reproducibility in Coagulation Assays, *AJCP*, **55**:561-564
- Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, *AJCP*, **59**:231-235

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BIBLIOGRAFIA

- Babson AL and Flanagan ML (1975) Quantitative One Stage Assays for Factors V and X, *AJCP*, **64**:817-819
- Hardisty RM, *et al.* (1962) A One Stage Factor VIII Assay and Its Use on Venous and Capillary Plasma, *Thrombosis et Diathesis Haemorrhagica*, 7:215-229
- Morse EE, *et al.* (1971) Automated Fibrinogen Determination, *AJCP*, **55**:671-676
- Elodi S, Katalin V, Hollan S (1978) Some Sources of Error in the One-Stage Assay of Factor VIII, *Haemostasis*, 7:1-9
- Thelin M (1968) Preparation and Standardization of a Stable AHF Plasm, *Thrombosis et Diathesis Haemorrhagica*, **19**:423
- Kirkwood TBL et al. (1977) Identification of Sources of Variation in Factor VIII Assay, *British Journal of Haematology*, **37**:559-568
- Goldenfarb MD (1971) Reproducibility in Coagulation Assays, *AJCP*, **55**:561-564
- Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, *AJCP*, **59**:231-235

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