# Thrombin Time

REF 5377  
REF 5392  
REF 5392H

Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 0SD, United Kingdom  
Tel: +44 (0)191 482 8440  
Fax: +44 (0)191 482 8442  
Email: info@helena-biosciences.com  
Web: www.helena-biosciences.com

HL-2-0444P 2016/01 (12)

## Thrombin Time

Instructions for use

en

### INTENDED PURPOSE

The Thrombin Time kit is intended for carrying out clot based haemostasis assays.

The Thrombin Time reagent contains bovine thrombin. The reagent can be used manually or on semi-automated and automated instruments. The test is commonly applied to detect various sources of interference with normal blood coagulation. Prolongation of the thrombin clotting time can be taken as a qualitative indication of abnormal fibrinogen levels (high or low), or the presence of interfering substances such as Fibrin Degradation Products (FDPs) or heparin. Quantitative evaluation of the possible causes of prolonged thrombin clotting time should be performed as follow-up studies, such as aPTT or chromogenic assay for heparin, Clauss fibrinogen, FDP determinations, heparin neutralisation by protamine sulphate or polybrene<sup>1</sup>, normal plasma mixing studies<sup>2</sup> or reptilase assay<sup>3</sup> to distinguish between hypofibrinogenaemia and FDP effects.

### WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

### COMPOSITION

Composition	Content	Description	Preparation
Thrombin Time	10 x 1 mL (REF 5392H)	Each vial contains a lyophilised preparation of bovine thrombin with buffers and stabilisers. The reconstituted reagent contains ≤10 NIH units/mL of thrombin. The reagent should be a white lyophilised plug.	Reconstitute the reagent with the recommended volume of purified water: 1 mL - REF 5392H 2 mL - REF 5392 5 mL - REF 5377
	10 x 2 mL (REF 5392)		Allow to stand for 5 minutes then mix gently by inversion and transfer to a plastic tube. Reconstituted reagent should be a clear, colourless solution.
	10 x 5 mL (REF 5377)		

Each kit contains Instructions For Use.

### ITEMS REQUIRED BUT NOT PROVIDED

Any high quality, mechanical, opto-mechanical or photo-optical coagulation instrument capable of performing thrombin clotting time testing may be used.

### STORAGE, SHELF-LIFE AND STABILITY

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Thrombin Time Reconstituted reagent should be stored at  $^{\circ}2$ – $^{\circ}8$ C and is stable for 14 days or at  $-20$ C for 1 month.  
Reconstituted reagent can also be stored on-board a Sysmex CA1500 photo-optical instrument for up to 7 days.

### SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at  $^{\circ}2$ – $^{\circ}8$ C or  $^{+18}$ – $^{+24}$ C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at  $-20$ C for 2 weeks or  $-70$ C for 6 months. Thaw quickly at  $^{\circ}37$ C prior to testing. Do not keep at  $^{\circ}37$ C for more than 5 minutes<sup>4</sup>.

### PROCEDURE

- Collect and prepare the blood sample according to the directions outlined in "Sample Collection And Preparation".
- Reconstitute the control plasma according to the package insert provided with each control.
- Prepare the reagent for use in the procedure according to the reconstitution instructions in the "Composition" section.
- Perform all tests in duplicate. Calculate the mean clotting time of the duplicate determinations to the nearest 0.1 seconds. Individual values should be within  $\pm 5\%$  of the mean value.

#### Manual Method

- Pipette 0.2 mL of the patient plasma or control plasma into a reaction tube.
- Incubate at  $^{\circ}37$ C for 3 minutes.
- Pipette 0.1 mL of Thrombin Time into the reaction tube containing patient or control plasma whilst simultaneously starting a timer.
- Record the time taken for clot formation to the nearest 0.1 second.

#### Automated Method

Refer to the appropriate instrument operator manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application notes.

### INTERPRETATION OF RESULTS

The results of the Thrombin Time test should be reported to the nearest 0.1 seconds. The normal range (usually the mean  $\pm 2$  standard deviations) should be established by each laboratory. Results outside the normal range should be considered abnormal and follow-up testing performed.

### LIMITATIONS

Expected values for the Thrombin Time test will vary from one laboratory to another depending on the techniques used. The method of clot detection, temperature, pH, collection technique, type of anticoagulant and time and method of sample storage are all very important. Plasma sample collection and storage conditions should be standardised and carefully controlled. Unexpected results should be confirmed by additional testing.

As well as the cause of elongated thrombin times indicated, a report has suggested that many systemic amyloidosis patients with bleeding complications may have a circulating inhibitor which prolongs the Thrombin Time<sup>5</sup>. Also, therapeutic levels of heparin may entirely abolish clotting in the Thrombin Time test, although neutralisation with protamine sulphate or polybrene should correct the Thrombin Time<sup>6</sup>.

### QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. Helena Biosciences Europe supplies the following controls available for use with this product:

REF 5186 Routine Control N  
REF 5187 Routine Control A  
REF 5183 Routine Control SA

### REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges.

### PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives:

#### Reproducibility

Sample	n	Intra-assay precision		Inter-assay precision	
		Clot formation (seconds)	CV (%)	Clot formation (seconds)	CV (%)
Normal	10	12.5	1.43	12.5	1.29

### BIBLIOGRAPHY

- Laposata *et al.* The Clinical Haemostasis Handbook, Yearbook Medical Publishers Inc., p219, 1989
- Thompson AR and Harker LA. Manual of Haemostasis and Thrombosis. 3rd Ed., F.A. Davis Co., p62, 1983
- DeMott WR. Laboratory Test Handbook, 2nd Ed., Jacobs D.S. et al Eds., Lexi-Comp Inc., p432-433, 1990
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Gastineau DA *et al.* (1991) Inhibitor of the Thrombin Time in Systemic Amyloidosis: A Common Coagulation Abnormality, *Blood*, 77:2637-40

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## BIBLIOGRAFIA

- Lapostata et al. The Clinical Haemostasis Handbook, Yearbook Medical Publishers Inc., p219, 1989
- Thompson AR and Harker LA. Manual of Haemostasis and Thrombosis. 3rd Ed., F.A. Davis Co., p62, 1983
- DeMott WR. Laboratory Test Handbook, 2nd Ed., Jacobs D.S. et al Eds., Lexi-Comp Inc., p432-433, 1990.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Gastineau DA et al. (1991) Inhibitor of the Thrombin Time in Systemic Amyloidosis: A Common Coagulation Abnormality, *Blood*, 77:2637-40

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