# Calibration Plasma

REF 5185

<span></span>	Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom
<span></span>	Tel: +44 (0)191 482 8440
<span></span>	Fax: +44 (0)191 482 8442
<span></span>	Email: <a href="mailto:info@helena-biosciences.com">info@helena-biosciences.com</a>
<span></span>	Web: <a href="http://www.helena-biosciences.com">www.helena-biosciences.com</a>

HL-2-0424P 2016/01 (14)

## Calibration Plasma Instructions for use

### INTENDED PURPOSE

The Calibration Plasma kit is intended for use as a calibration material.

Calibration Plasma may be used as a reference plasma when assaying for Factors II, V, VII, VIII, IX, X<sup>1</sup>, XI, XII<sup>2</sup>, fibrinogène<sup>3</sup>, von Willebrand Factor, antigenic and functional Protein C and Protein S (total and free), as well as the chromogenic assays including Antithrombin (AT) III, Protein C, Factor VIII and Plasminogen. Calibration Plasma is used in factor assays and other testing in the same manner as a fresh plasma pool. Factor II, VII, VIII, IX and X values and the chromogenic Factor VIII, AT III and Protein C values are traceable to World Health Organisation standards to ensure the utmost credibility in values stated<sup>4</sup>. The reference plasma should be used to gauge internal factors in each laboratory’s system. The reference plasma should not be used to determine normal ranges since normals vary from population to population.

### WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

Blood products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of:
Hepatitis B Antigen (HbsAg)
HIV 1 antibody
HIV 2 antibody
HCV antibody
However they should be handled with the same precautions as a human patient sample.

### COMPOSITION

Component	Content	Description	Preparation
Calibration Plasma	10 x 1 mL	Calibration Plasma is prepared from a frozen pool of citrated plasma from healthy donors and is buffered and lyophilised to ensure stability of all plasma constituents <sup>5</sup> .	Reconstitute with 1 mL of distilled or deionised water. Swirl gently. Allow product 20 minutes for complete dissolution.

Each kit contains instructions for use.

Each kit contains a lot specific reference values insert.

### ITEMS REQUIRED BUT NOT PROVIDED

Calibration Plasma may be used when performing tests on any mechanical or photo-optical coagulation instrument in conjunction with all suitable, commercial reagents.

### STORAGE, SHELF-LIFE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label. Values for Factor VIII, von Willebrand factor and ristocetin co-factor are stable for 2 hours at \*2 –8°C. All other factors are stable for 4 hours at \*2 –8°C. The unreconstituted Calibration Plasma must appear as a light yellow, dry plug. If any unusual conditions are noted, the customer should notify Helena Biosciences Europe, before using.

### SAMPLE COLLECTION AND PREPARATION

Not applicable.

### PROCEDURE

Refer to the appropriate instrument Operator Manual for detailed instructions or contact Helena Biosciences Europe for instrument-specific application guides.

### INTERPRETATION OF RESULTS

The assayed percent activities of the various coagulation factors should be taken from the "Reference Value" column when Calibration Plasma is used to determine the laboratory’s standard curves. Ensure the lot number printed on this assay sheet is the same as that on the vial of Calibration Plasma to be used.

### LIMITATIONS

The results obtained with Calibration Plasma depend on several factors strongly associated with instrumentation, types of reagents, deficient substrates and laboratory to laboratory variations<sup>6,7,8</sup>. Each laboratory should establish an expected range for the particular instrument-reagent system.

### QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Helena Biosciences Europe supplies the following controls available for use with this product:

REF 5301 Speciality Assayed Control N

REF 5302 Speciality Assayed Control A

### REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range.

### PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena Biosciences Europe or their representatives using an opto-mechanical coagulation instrument. Each laboratory should establish its own performance data.

<b>Reproducibility</b>			
<b>Inter-batch precision</b>			
<i>Parameter</i>	<i>n</i>	<i>Mean</i>	<i>CV (%)</i>
Fibrinogen (g/L)	5	3.0	3.6
Facteur IX (%)	5	124.7	1.9
Protéine S (%)	5	98.5	2.8

### BIBLIOGRAPHY

- Babson AL and Flanagan ML (1975) Quantitative One Stage Assays for Factors V and X, *AJCP*, **64**: 817-819.
- Hardisty RM *et al.* (1962) A One Stage Factor VIII Assay and Its Use on Venous and Capillary Plasma. *Thrombosis et Diathesis Haemorrhagica*, **7**:215-229.
- Morse EE *et al.* (1971) Automated Fibrinogen Determination, *AJCP*, **55**:671-676.
- Elodi S *et al* (1978) Some Sources of Error in the One-Stage Assay of Factor VIII, *Haemostasis*, **7**:1-9.
- Thelin M (1968) Preparation and Standardization of a Stable AHF Plasma, *Thrombosis et Diathesis Haemorrhagica*, **19**:423.
- Kirkwood TBL *et al.* (1977) Identification of Sources of Variation in Factor VIII Assay, *British Journal of Haematology*, **37**:559-568.
- Goldenfarb MD (1971) Reproducibility in Coagulation Assays, *AJCP*, **55**:561-564.
- Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, *AJCP*, **59**:231-235.

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