REF 5374

REF 5378

 	
	Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom
	Tel: +44 (0)191 482 8440
	Fax: +44 (0)191 482 8442
	Email: info@helena-biosciences.com
	Web: www.helena-biosciences.com

HL-2-0651P 2015/10 (10)

Claus Fibrinogen (Thrombin only)	en
Instructions for use	

INTENDED PURPOSE

The Claus Fibrinogen (Thrombin only) kit is intended for carrying out clot based haemostasis assays.

Claus¹ developed a simple method for the quantitative determination of fibrinogen by measuring the clotting time of dilute plasma after the addition of Thrombin (>30 NIH units/mL). This clot time is proportional to the fibrinogen concentration. Levels of fibrinogen can increase as a result of inflammation, pregnancy or oral contraceptive use². Decreased levels can be found in certain states such as liver disease and DIC. Congenital deficiencies include afibrinogenaemia (no detectable fibrinogen), hypofibrinogenaemia (<1 mg/mL) and dysfibrinogenaemia (abnormal fibrinogen molecule). Claus Fibrinogen (Thrombin only) is intended for use in the quantitative determination of fibrinogen in human plasma using the Claus¹ Method.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

COMPOSITION

Component	Content	Description	Preparation
Thrombin: 100 NIH/mL	10 x 2 mL (REF 5374) <p>10 x 5 mL (REF 5378)</p>	Each vial contains approximately 100 NIH units/mL bovine thrombin with stabilisers.	Reconstitute each vial with: 2 mL (REF 5374) 5 mL (REF 5378) of purified water. Take care when pipetting to avoid contamination. Swirl gently and allow to stand for 15 minutes. Mix well immediately before use. Do not shake.

Each kit contains Instructions For Use.

ITEMS REQUIRED BUT NOT PROVIDED

REF 5379	Fibrinogen Calibrator
REF 5185	Calibration Plasma
REF 5375	Owren's Buffer

Helena Biosciences Europe also provides the complete fibrinogen assay kit, REF 5376, Claus Fibrinogen 100.

STORAGE, SHELF-LIFE AND STABILITY

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

Thrombin: Once reconstituted, the reagent is stable for 8 hours at +15 –+30°C, 1 week at *2 –+8°C or 1 month at -20°C.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at *2 –+8°C or *18 –+24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at *37°C prior to testing. Do not keep at *37°C for more than 5 minutes⁵.

PROCEDURE

Manual Method

Prepare all reagents as instructed under "composition".

- Standard Curve Preparation:
 - A new standard curve should be used when reagent lot numbers change or if expected values in quality control testing change. Prepare the following dilutions in Owren's Buffer (mix without shaking):

Tube	Dilution	Fibrinogen Calibrator (mL) (REF 5379 / REF 5185)	Buffer (mL) (REF 5375)
1	1 + 4	0,2	0,8
2	1 + 9	0,1	0,9
3	1 + 19	0,1	1,9
4	1 + 29	0,1	2,9
5	1 + 39	0,1	3,9

- Patient Sample Preparation:
 - Prepare 1+9 dilutions of the patient plasma or control plasma (REF 5185) in Owren's Buffer.
 - Mix without shaking.
- Testing:
 - Perform all tests in duplicate:
 - Pipette 0,2 mL of standard, patient or control dilution into a reaction tube and incubate at +37°C for 2 minutes.
 - Add 0,1 mL of thrombin (*15 –+30°C).
 - Determine the clot time to the nearest 0.1 seconds.
 - Plot mean Standard Clot Time (Y-axis) versus Fibrinogen Level (X-axis) on the Fibrinogen graph paper. A straight line should be obtained. Assign the Fibrinogen Calibrator reference value to the 1+9 dilution to allow direct interpolation of patient and control values from the standard curve.

Automated Method

Refer to the appropriate instrument Operator Manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application guides.

INTERPRETATION OF RESULTS

Expected values for fibrinogen in healthy adults are 150-350 mg/dL (1.5-3.5 g/L)^{4,5}.

LIMITATIONS

Heparin levels >0.6 units/mL and fibrinolytic degradation products >100 mg/mL may cause falsely low fibrinogen quantitation. If values fall outside the standard curve values for the patient samples, re-assay using an appropriate dilution to bring values into the standard range.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. Helena Biosciences Europe supply the following controls available for use with this product:

REF 5301	Speciality Assayed Control N
REF 5302	Speciality Assayed Control A
REF 5186	Routine Control N
REF 5187	Routine Control A
REF 5183	Routine Control SA

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range.

PERFORMANCE CHARACTERISTICS

Helena Biosciences Europe or their representatives have determined the following performance characteristics as a guideline. Each laboratory should establish its own performance data. The manual Claus Fibrinogen is designed to give a linear calibration from 1.5 - 6.5 g/L.

Reproducibility					
Intra-assay precision	Inter-assay precision				
<i>Fibrinogen (g/L)</i>	<i>n</i>	<i>CV (%)</i>	<i>Fibrinogen (g/L)</i>	<i>n</i>	<i>CV (%)</i>
3.64	10	4.87	3.61	10	3.78

BIBLIOGRAPHY

- Claus A (1957) Gerinnungsphysiologische Schnell-methode zur bestimmung des Fibrinogens, *Acta Haematol*, **17**:237-246
- Shaw TS (1977) Assays for Fibrinogen and its Derivatives CRC, *Crit. Rev. Clin. Lab. Sci*, **8**:145-192.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Scully RE *et al.* (1980) Normal Reference Laboratory Values, *New England Journal Medicine*, **302**(37):37-48.
- Okuno T, Selenko V (1972) Plasma fibrinogen determination by automated thrombin time, *American Journal of Medical Technology*, **38**(6):196-201

Claus Fibrinogen (Thrombin only)	fr
Fiche technique	

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BIBLIOGRAFIA

- Clauss A (1957) Gerinnungsphysiologische Schnell-methode zur bestimmung des Fibrinogens, *Acta Haematol*, **17**:237-246
- Shaw TS (1977) Assays for Fibrinogen and its Derivatives CRC, *Crit. Rev. Clin. Lab. Sci*, **8**:145-192.
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Scully RE *et al.* (1980) Normal Reference Laboratory Values, *New England Journal Medicine*, **302**(37):37-48.
- Okuno T, Selenko V (1972) Plasma fibrinogen determination by automated thrombin time, *American Journal of Medical Technology*, **38**(6):196-201

Clauss Fibrinogen (Thrombin only) Instrucciones de uso	es
---------------------------------------------------------------	-----------

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