REF 5485

 	
	Helena Biosciences Europe, Queensway South, Team Valley Trading Estate, Gateshead, Tyne and Wear, NE11 OSD, United Kingdom
	Tel: +44 (0)191 482 8440
	Fax: +44 (0)191 482 8442
	Email: info@helena-biosciences.com
	Web: www.helena-biosciences.com

HL-2-0722P 2016/01 (8)

DRVVT Confirm Instructions for use	en
INTENDED PURPOSE	
The DRVVT Confirm kit is intended for carrying out clot based haemostasis assays.	

Lupus Anticoagulants (LA) are antibodies of the IgG and IgM type which are directed against a variety of anionic phospholipids. The presence of LA in plasma is increasingly associated with a variety of haemostatic problems such as recurrent fetal loss, thrombocytopenia, unexplained thrombosis and neurological disorders. LA prolongs phospholipid dependant *in vitro* clotting assays such as the activated partial thromboplastin time (aPTT)

The Helena Biosciences Europe DRVVT Confirm kit is intended for the qualitative confirmation of LA in human plasma. The reagent is designed to be used in conjunction with the DRVVT Screen Kit (REF 5484) to discriminate between LA, factor deficiencies (II, V or X) or other inhibitors. If the clot time of the patient samples with the DRVVT Screen procedure are greater than 3 standard deviations above the mean of the normal range and are not corrected by mixing studies, a lupus anticoagulant may be present. Under these circumstances, samples should be re-tested using the DRVVT Confirm Reagent.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in vitro* diagnostic use only – DO NOT INGEST. Wear appropriate personal protective equipment when handling all kit components. Refer to the product safety declaration for the link to appropriate hazard and precautionary statements where applicable. Dispose of components in accordance with local regulations.

Component	Content	Description	Preparation
DRVVT Confirm	10 x 1 mL	Each vial contains a proprietary mixture of Russell’s Viper Venom co-lyophilised with calcium chloride and concentrated phospholipid.	Reconstitute each vial with 1 mL of distilled / deionised water. Allow to stand for 10 minutes and mix well before use. Do not shake.

Each kit contains Instructions For Use.

ITEMS REQUIRED BUT NOT PROVIDED

Optional Items:

REF 5484 DRVVT Screen

STORAGE, SHELF-LIFE AND STABILITY

Unopened reagents are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

DRVVT Reconstituted vials are stable for 24 hours at *15 –*30°C, 5 days at *2 –*8°C or 2 weeks at -20°C. The reagent should be frozen in plastic test tubes and thawed at *37°C before use.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at *2 –*8°C or *18 –*24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at *37°C prior to testing. Do not keep at *37°C for more than 5 minutes¹. If freezing, double centrifugation of the sample is recommended to ensure that the sample is platelet poor. Transfer the plasma following the initial centrifugation to a non-activating plastic tube using a plastic pipette, then re-centrifuge the plasma for an additional 10 minutes at a higher speed (>2500 x g). When aliquoting to a secondary tube, take care to not include the residual platelets that may have collected at the bottom of the centrifuge tube².

PROCEDURE

Manual Method

- Pre-warm sufficient reconstituted reagent to *37°C.
- Pipeter 0.2 mL of patient or control plasma into a reaction tube. Incubate at *37°C for 2 minutes.
- Add 0.2 mL of pre-warmed DRVVT Confirm reagent and start a timer.
- Measure the clot formation time to the nearest 0.1 seconds.
- Calculate the normalised ‘DRVVT Confirm’ ratio as:

Patient DRVVT Confirm Clot Time / Mean Normal DRVVT Confirm Clot Time

Automated Method

Refer to the appropriate instrument operator manual for detailed instructions or contact Helena Biosciences Europe for instrument specific application guides.

INTERPRETATION OF RESULTS

Results are best expressed as a normalised ratio relative to the mean normal clot time obtained by each laboratory. It is recommended that like for like sample types are used when calculating a normalised ratio. Both DRVVT Screen and DRVVT Confirm results can be ‘normalised’ in this way, reducing the effect of instrument variability and potentially improving discrimination between weak positive LA and normal samples. Results of the mixing tests can be treated in the same way.

If patient samples with suspected LA (by DRVVT Screen testing) have a DRVVT Confirm clot time less than 3 SDs from the mean normal clot time for this kit, LA is strongly indicated. Mixing studies would also give normal clot times in this situation. Samples with an abnormal DRVVT Screen clot time and abnormal DRVVT Confirm clot time which are both normalised by mixing studies are most likely factor deficient (II, V or X). Samples with abnormal DRVVT Screen clot time and abnormal DRVVT Confirm clot time in which only the DRVVT Confirm clot time is normalised in mixing studies suggests LA plus factor deficiency. Samples with abnormal DRVVT Screen clot time and abnormal DRVVT Confirm clot time which are not corrected by mixing studies indicate other inhibitors. The results of DRVVT Screen and DRVVT Confirm testing when expressed as the ‘normalised’ ratio can also be used to indicate the level of LA present:

(DRVVT Screen Ratio)/(DRVVT Confirm Ratio) > 2.0 Strong LA
(DRVVT Screen Ratio)/(DRVVT Confirm Ratio) 1.5 - 2.0 Moderate LA
(DRVVT Screen Ratio)/(DRVVT Confirm Ratio) 1.2 - 1.5 Weak LA

LIMITATIONS

Plasma deficiencies of Factors II, V or X may lead to abnormal results in neat plasma. Mixing studies should correct this. Plasma from patients with the following may give abnormal results when the plasma is tested neat, and these samples may not correct in mixing studies: heparin (>1 U/mL), oral anticoagulants, disseminated intravascular coagulation (DIC). Care must be taken to remove residual platelets from plasma by filtration or centrifugation, as platelet derived phospholipid can interfere with the test.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid. Helena Biosciences Europe supplies the following controls available for use with this product:

REF 5486 LA Positive Control S
REF 5186 Routine Control N

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own reference ranges. The normal reference range (mean ± 3SDs) determined at Helena Biosciences Europe for the DRVVT Confirm test was 32.4 ± 6.0 seconds (range 26.4-38.4 seconds).

PERFORMANCE CHARACTERISTICS

Each laboratory should establish its own performance data. Within run and between run CVs are expected to be <5%.

BIBLIOGRAPHY

- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5
- Pengo *V et al* (2009) Update of the guidelines for lupus anticoagulant detection, *J Thromb Haemost*, **7**: 1737-40

DRVVT Confirm Fiche technique	fr
---	-----------

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